



**COVALENT**  
**METROLOGY**

**Welcome**

# LASER ABLATION INDUCTIVELY COUPLED PLASMA MASS SPECTROSCOPY: NOT JUST ROCKS

SPEAKER:

Lucas Smith

Director of Business Development  
for the Americas,  
Teledyne CETAC

February 10, 2022 | 11am PT



## COVALENT ACADEMY


Advancements in  
Instrumentation Series

Episode 30



# COVALENT METROLOGY

Silicon Valley-based analytical labs and platform delivering quality data and expert analysis for advanced materials and device innovation



## Comprehensive Solutions Stack

*50+ cutting-edge instruments, offering 100+ Techniques*

*Analytical Services*

*Advanced Modeling*

*Method Development*

*Temp. Staffing Solutions*



## Affordable and Fast

*Fast Turnaround Times, No Expedite Fees*

*Volume Savings*

*Instant Access to Data and Reports in Secure Portal*



## Flexible Business Model

*Custom Consulting Solutions and Certified Onsite Support*

*Training and Certification on Instrumentation*

*Co-op and Tool-Share Opportunities*

*Laboratory Audits*



## Rich Network of Partnerships

*Partner to World's Leading Instrument Manufacturers and Labs*

*Expanding Instrumentation, Lab Connections and Learning*



## Who We Are, Who We Serve

*500+ Clients, 40-60 Added / Quarter*

*50+ People, 14 PhDs*

*Cutting-edge Analytical Capabilities*

*Lab Location: Sunnyvale, CA*

# Covalent Technical Groups and Organization

4

## PCBA, Semiconductor, and Electronic Device Metrology & Failure Analysis

- DPA / Mechanical Cross-section
- Dye & Pry Test
- EBIC / OBIC failure analysis
- Hot Spot Detection
- IR Imaging / Emission Microscopy
- NIR Imaging
- Root-Cause Failure Analysis

## Electron Microscopy and Scanning Probe Microscopy

- AFM & Advanced AFM Modes (EFM, KPFM, MFM, PFM)
- Scanning Acoustic Microscopy (SAM)
- SEM (+ EDS)
- FIB-SEM ( + EDS)
- S/TEM ( + EDS / + EELS )
- Nano-indent / Nano-scratch

## Optical Microscopy & Spectroscopy

- Chromatic Aberration
- Digital Optical Microscopy
- FTIR and ATR-FTIR
- Laser Scanning Confocal Microscopy
- Spectral Ellipsometry
- UV-Vis-NIR Spectroscopy
- White Light Interferometry

## X-Ray Characterization

- X-Ray Diffraction (XRD)
- X-Ray Reflectometry (XRR)
- Micron-spot ED-XRF
- WDXRF
- Micro-computed X-ray Tomography (Micro-CT)
- 2D X-ray Inspection & X-ray Radiography

## Elemental / Chemical Composition Analysis

- EPMA
- GD-OES
- GC-MS
- ICP-MS and LA-ICP-MS
- Raman Microscopy & Spectroscopy
- NMR (1D or 2D; solid / liquid)

## Particle Analysis

- Dynamic Light Scattering (DLS)
- Laser Diffraction Particle Size Analysis (PSA)
- Particle Zeta Potential

## Material Property Characterization

- DSC
- DMA & TMA
- Rheometry
- TGA
- Surface Zeta Potential
- Porometry / Porosity
- Gas Adsorption
- Gas Pycnometry
- Foam Density
- Tap Density

## Surface Spectroscopy Analysis

- Dynamic-SIMS
- ToF-SIMS (Static-SIMS)
- Ion Scattering Spectroscopy (ISS)
- Ultraviolet Photoelectron Spectroscopy (UPS)
- X-ray Photoelectron Spectroscopy (XPS)



- Teledyne CETAC is a global leader in the development of advanced sample-introduction and sample-handling equipment for elemental analysis
- Partnership with Covalent announced in Fall, 2021 to:
  - Jointly investigate and develop new applications for Teledyne's element analysis tools
  - Bring cutting-edge analytical chemistry solutions to Covalent clients
- First Teledyne instrument installed at Covalent is the *Iridia Laser Ablation System*

## Other Covalent Partners



## Lucas Smith

Director for Business Development for the Americas,  
Teledyne CETAC





---

- Degrees in Microbiology and Biochemistry from The University of Michigan (Go Blue!)
- Manages laser opportunities and customer relationships throughout the Western Hemisphere
- Member of the Teledyne team for just over 7 years
- Prior to joining Teledyne Lucas held a variety of roles at
  - Pfizer Pharmaceuticals
  - LECO Corporation
- Performed mobile analytics for the US Navy & Airforce




# Content

- What is laser ablation ICP-MS?
- Why laser ablation (advantages vs disadvantages)
- Laser ablation vs other analytical techniques
- Spatial Analysis
- Volumetric Ablation
- Application examples
- Q&A




# What is Laser Ablation?

- A laser ablates a solid sample and delivers material directly to the ICPMS for elemental analysis
- Eliminates dissolution step and most sample prep
- Maintain spatial resolution




# Why Laser Ablation?

- Rapid, sensitive multi-element detection at sub-ppm levels.
- Eliminates (most) labor-intensive sample preparation
- Can sample virtually any material (including liquids)
- Flexible method development
- High spatial resolution of microfeatures
- Sample mapping and depth profiling capabilities
- No contamination of reagents from digestion/dilutions steps




# Challenges of the Technique


- Historically difficult acquiring standards
- Upper limit in sensitivity due to detector saturation
- (micro)Destructive technique
- Initial cost




# Laser Ablation VS Other Techniques



# What is Laser Ablation ICP-MS?





# Spatial analysis




# Spatial analysis

- OD – Single Spot - this is (still..!) the typical LA-ICP-MS analysis
- Suitable for bulk analysis, depth profiling, isotope ratio, etc.



# Spatial analysis




Line

- 1D – Transect
- Useful for the study of diffusion profiles, bulk analysis of homogeneous samples (e.g., Particles trapped on filters, etc.)


Ogrizek, M., Jaćimović, R., Šala, M., & Kroflič, A. (2021). No more waste at the elemental analysis of airborne particulate matter on quartz fibre filters. *Talanta*, 226.



# Spatial analysis




- 2D – the spatial distribution of elements (isotopes) in samples



Metamorphic garnet from the Eastern Carpathians (Romania)  
Courtesy of Dr. Gavril Sabau, Geological Institute of Romania

# Spatial analysis



- 3D – the spatial distribution of elements (isotopes) in 3D



Zircon from the Ditrau alkaline massif (Eastern Carpathians, Romania).  
Courtesy of Dr. Gavril Sabau, Geological Institute of Romania


# Data Collected at 375Hz

## Mouse eye

- 3.0mm diameter
- Unpublished data

**10  $\mu\text{m}$  resolution  
375 Hz**

- Agilent 8900 ICP-MS
- Cobalt *Long Pulse* cell
- Parameters optimized and images generated with Teledyne HDIP software
- **5 elements measured**
  - $\text{C}^{13}$ ,  $\text{Fe}^{57}$ ,  $\text{Cu}^{63}$ ,  $\text{Zn}^{66}$ ,  $\text{Pb}^{208}$




50 minutes

**N** NATURAL HISTORY MUSEUM


Sheffield Hallam University

# Horizontal & Vertical Data Transects



Transects extracted with Teledyne HDIP software

# Volumetric Ablation – for Depth Profiling and Thin Film Analysis




FC-1 zircon

3  $\text{Jcm}^{-2}$  (measured at the sample surface)

10 Hz

20  $\mu\text{m}$  spot

1, 5, 10, 20, 50, 100 and 157 pulses




FC-1 zircon

3  $\text{Jcm}^{-2}$

10 Hz

20  $\mu\text{m}$  spot


1 pulse




# Volumetric Ablation – for Depth Profiling and Thin Film Analysis

Shots	Depth	Volume	Expected
1	0.04	13.26	13.26
5	0.25	91.56	66.28
10	0.5	171.0	132.6
20	1	349.0	265.1
50	2.7	708.5	662.8
100	5.55	1070	1326
157	8.1	1954	2081

Zircon FC-1, 20 μm spot, 3 J/cm<sup>2</sup>






Zircon FC-1, 20 μm spot, 3 J/cm<sup>2</sup>



# Thin Film Analysis

Ion signal evolution for 5 individual laser pulses


Sample: 6 nm Nd internal layer




High energy laser pulses can remove thin layers within a few laser pulses




# Elemental Imaging – Basics of Map/Image Generation



**Line-by-line ablation**




**Transient ion traces**



**Ion intensity image**

# Elemental Imaging





	A	B	C	D	E	F	G	H	I	J	K	L	M
8	0.096912	0.114887	0.583607	0.660054	0.852923	0.218161	0.397852	0.78179	0.633538	0.383446	0.096522	0.502717	1.18885
9	0.146324	0.161749	0.464136	0.69083	0.623064	0.204618	0.404379	0.787971	0.633808	0.384757	0.198939	0.33847	1.19359
10	0.16358	0.155945	0.601307	0.395216	0.692919	0.485034	0.198672	0.556445	0.57554	0.364166	0.13546	0.333279	1.07667
11	0.235245	0.048706	0.774869	0.412034	1.03772	0.49413	0.396634	0.292094	0.385738	0.438952	0.193673	0.304141	1.02072
12	0.39712	0.14578	0.842969	0.40835	0.946479	0.643883	0.193839	0.556434	0.453219	0.189042	0.234477	0.263653	0.662822
13	0.228604	0.079104	1.19217	0.830759	1.16044	0.536759	0.131089	0.130565	0.244766	0.250621	0.202918	0.497544	0.532865
14	0.259311	0.237749	0.864093	0.823291	1.2092	0.696092	0.162978	0.415797	0.220872	0.089312	0.210114	0.775358	0.56165
15	0.390668	0.190997	1.06221	0.999354	1.30664	0.866227	0.227986	0.176978	0.086577	0.13053	0.248586	0.472006	0.35736
16	0.104634	0.409608	1.61818	1.43177	1.31183	0.545197	0.293767	0.264986	0.16443	0.072054	0.091039	0.754845	0.253603
17	1.16384	0.385957	1.65124	0.798049	1.56494	0.66774	0.415711	0.573101	0.046124	0.074338	0.611663	0.610238	0.226235
18	0.369656	0.715347	1.53221	6.00476	1.38987	1.05974	0.666762	1.13793	0.026259	0.211942	0.457764	0.350029	0.116701
19	0.352381	1.07845	4.65948	2.63666	1.17151	1.08914	0.720548	1.27075	0.068776	0.145108	0.386675	0.258985	0.04538
20	0.505702	1.10533	1.41073	1.47368	1.21054	0.930509	0.639065	0.895207	0.026999	0.050367	0.28623	0.430147	0.12626
21	0.599119	0.896266	0.525086	1.20172	1.30803	0.803165	0.563338	0.772326	0.156519	0.096277	0.404549	0.324726	0.103211
22	0.414209	0.801911	0.843391	1.09828	1.25331	0.985016	0.449396	0.53338	0.099202	0.072237	0.305811	0.368388	0.133242
23	0.366787	0.768113	0.531492	1.07642	1.11176	0.987023	0.758482	0.611338	0.056195	0.045279	0.151004	0.23329	0.064937
24	0.458827	0.845649	0.824528	1.03083	1.07516	1.18947	0.514221	0.293982	0.159779	0.113608	0.22589	0.171253	0.041933
25	0.716285	0.663331	0.689056	0.820868	1.02739	1.16858	0.718875	0.474105	0.100959	0.083564	0.164511	0.105662	0.083173
26	0.819503	1.24763	1.42895	0.707344	1.52214	1.09088	0.763535	1.13511	0.033873	0.142139	0.167168	0.099386	0.043056
27	0.912334	2.01484	1.20841	0.724997	1.44369	1.09532	0.917739	1.39346	0.048793	0.14082	0.097433	0.074208	0.03787
28	1.15607	3.04507	1.26789	0.62716	1.25552	0.793643	1.00461	0.989637	0.03223	0.085011	0.113863	0.018183	0.098037
29	1.15575	2.83549	1.06188	0.610768	1.54787	0.909811	0.55543	1.672	0.064739	0.10096	0.171191	0.05557	0.056548
30	1.12294	1.72703	0.717315	0.934489	1.30029	1.02686	0.785505	1.44898	0.024459	0.053282	0.110676	0.115487	0.096095
31	1.67334	1.15917	0.571676	0.631869	1.5819	0.646488	0.656882	1.15551	0.027585	0.079305	0.118181	0.086831	0.154575
32	1.6434	0.868305	0.285532	0.760374	1.91634	0.746576	0.729745	1.14162	0.033734	0.076443	0.050522	0.097287	0.204627
33	1.65086	0.671309	0.166098	0.58059	1.78824	0.983115	0.774283	1.11937	0.034717	0.162403	0.054783	0.087952	0.131606
34	1.30019	0.857228	0.160557	0.581421	1.67726	2.02466	1.25894	0.751958	0.016627	0.103778	0.056284	0.066515	0.162199
35	1.00326	0.564484	0.228819	0.760232	1.27207	1.84559	1.43669	0.496919	0.041127	0.173992	0.057838	0.080648	0.181131
36	0.775591	0.41202	0.417392	0.642002	0.946925	1.4467	1.69969	0.354116	0.054937	0.128461	0.08167	0.083234	0.138208
37	0.662539	0.370612	0.532953	0.636069	0.983973	1.55404	1.32917	0.202616	0.021751	0.119431	0.072847	0.085618	0.138608
38	0.463375	0.31527	0.358647	0.586981	1.13203	1.17624	1.12798	0.188072	0.03207	0.065468	0.060936	0.055543	0.088213
39	0.540517	0.30932	0.42259	0.42403	1.23246	0.904797	0.651712	0.114582	0.110343	0.144098	0.092224	0.05494	0.060053
40	0.38211	0.329588	0.45556	0.863497	0.902195	0.664446	0.376653	0.074541	0.080114	0.064662	0.094256	0.129623	0.165235
41	0.421603	0.392585	0.279268	0.932795	0.900197	0.355988	0.44882	0.091691	0.170302	0.099841	0.138096	0.078876	0.125028
42	0.390633	0.381289	0.360255	0.811816	0.628627	0.361066	0.198197	0.126083	0.085758	0.11067	0.039471	0.077577	0.05944
43	0.347463	0.272831	0.298817	0.742719	0.252145	0.301427	0.175459	0.117354	0.067837	0.076817	0.032003	0.072841	0.083575



Courtesy of Dr. Gavril Sabau, Geological Institute of Romania






# Elemental Imaging

- Red blood cells of a patient after cisplatin chemotherapy (Courtesy of Dr. Sarah Theiner, University of Vienna)
- 2  $\mu\text{m}$  resolution
- 100 Hz
- Dosage: 4
- Imaging Time: 20 min



# Elemental Imaging


- Mouse Embryo (unpublished data)
- 5  $\mu\text{m}$  thickness, deparaffinized prior to analysis
- 8 hours analysis time for 18 samples
- 19 $\mu\text{m}$  Spot
- Dosage: 1
- Energy Density: 0.93 J  $\text{cm}^2$
- Scan Speed: 2 mm/s
- Rep Rate: 100Hz





# Elemental Imaging

- Barley root section
- 1  $\mu\text{m}$  resolution,
- 300 Hz
- 1 pulse/pixel
- Total Imaging Time: 3.2 min

Root4\_Mg\_1 $\mu\text{m}$



# HDIP Software – Powerful Simplicity (3D Imaging)



Three-Dimensional Reconstruction of the Tissue-Specific Multi-elemental Distribution within *Ceriodaphnia dubia* via Multimodal Registration Using Laser Ablation ICP-Mass Spectrometry and X-ray Spectroscopic Techniques. Stijn J. M. Van Malderen et Al, DOI: [10.1021/acs.analchem.7b00111](https://doi.org/10.1021/acs.analchem.7b00111)

# Data Collected at 400Hz




- Wheat grain (unpublished data)
- Mounted in epoxy resin and polished to give a flat surface
- Interested in zinc uptake
  - Concentration & location

# Data Collected at 400Hz


Carbon<sup>13</sup>

<sup>13</sup>C<sup>+</sup>



Zinc<sup>66</sup>

<sup>66</sup>Zn<sup>+</sup>



50 minutes

- **2 elements measured:** Carbon and Zinc
- **15 µm** circular spot; **400 Hz**
- Cobalt *Long Pulse* cell
- Agilent 8900 ICP-MS
- Parameters optimized and images generated with Teledyne HDIP software




**HIDIP 1.6**

Mass Spectrometry Data Analysis

Designed by scientists for scientists. Available for Teledyne Photon Machines laser systems.

# Elemental Imaging

- Buckwheat section (unpublished data).  
Courtesy of Dr. Martin Sala, National  
Institute of Chemistry Slovenia
- 10  $\mu\text{m}$  resolution
- 300 Hz
- 4 pulse/pixel







# Bulk Analysis

	Test 1 (10 Hz; Total Flow: 0.65 LPM)				Test 2 (50 Hz; Total Flow: 0.75 LPM)				Test 3 (50 Hz; Total Flow: 1.05 LPM)			
	[7Li]+	[89Y]+	[139La]+	[238U]+	[7Li]+	[89Y]+	[139La]+	[238U]+	[7Li]+	[89Y]+	[139La]+	[238U]+
Replicate 1 mean [cts]	2.62E+03	7.03E+03	8.89E+03	8.66E+03	8.81E+03	2.65E+04	3.20E+04	2.75E+04	8.81E+03	2.71E+04	3.25E+04	2.90E+04
Replicate 2 mean [cts]	2.66E+03	7.00E+03	8.74E+03	8.43E+03	8.97E+03	2.65E+04	3.22E+04	2.88E+04	8.84E+03	2.68E+04	3.22E+04	2.88E+04
Replicate 3 mean [cts]	2.71E+03	6.95E+03	8.68E+03	8.70E+03	8.98E+03	2.65E+04	3.24E+04	2.85E+04	9.09E+03	2.67E+04	3.24E+04	2.86E+04
Replicate 4 mean [cts]	2.74E+03	7.02E+03	8.69E+03	8.53E+03	9.13E+03	2.67E+04	3.22E+04	2.92E+04	9.03E+03	2.69E+04	3.25E+04	2.88E+04
Replicate 5 mean [cts]	2.73E+03	6.95E+03	8.85E+03	8.53E+03	9.16E+03	2.67E+04	3.26E+04	2.88E+04	9.06E+03	2.66E+04	3.22E+04	2.87E+04
Replicate 6 mean [cts]	2.81E+03	6.93E+03	8.62E+03	8.56E+03	9.09E+03	2.65E+04	3.24E+04	2.91E+04	9.06E+03	2.68E+04	3.24E+04	2.85E+04
Replicate 7 mean [cts]	2.72E+03	7.00E+03	8.81E+03	8.54E+03	9.15E+03	2.67E+04	3.25E+04	2.85E+04	8.90E+03	2.68E+04	3.25E+04	2.87E+04
Replicate 8 mean [cts]	2.68E+03	6.99E+03	8.78E+03	8.56E+03	8.96E+03	2.66E+04	3.22E+04	2.81E+04	8.95E+03	2.68E+04	3.23E+04	2.87E+04
Replicate 9 mean [cts]	2.65E+03	6.94E+03	8.75E+03	8.48E+03	8.98E+03	2.65E+04	3.22E+04	2.86E+04	8.92E+03	2.68E+04	3.25E+04	2.89E+04
Replicate 10 mean [cts]	2.65E+03	7.10E+03	8.77E+03	8.65E+03	8.92E+03	2.66E+04	3.22E+04	2.83E+04	8.95E+03	2.71E+04	3.26E+04	2.90E+04
Mean signal [cts]	2.70E+03	6.99E+03	8.76E+03	8.56E+03	9.02E+03	2.66E+04	3.23E+04	2.85E+04	8.96E+03	2.68E+04	3.24E+04	2.88E+04
Median signal [cts]	2.58E+03	6.83E+03	8.49E+03	8.21E+03	8.98E+03	2.65E+04	3.22E+04	2.85E+04	8.93E+03	2.68E+04	3.23E+04	2.86E+04
Measurement uncertainty mean (1SD) [cts]	1.65E+01	1.59E+01	2.42E+01	2.54E+01	3.41E+01	2.61E+01	5.31E+01	1.46E+02	2.89E+01	5.24E+01	3.85E+01	5.18E+01
Measurement uncertainty replicates (1SD) [cts]	5.22E+01	5.01E+01	7.64E+01	8.02E+01	1.08E+02	8.26E+01	1.68E+02	4.62E+02	9.12E+01	1.66E+02	1.22E+02	1.64E+02
Absolute total combined uncertainty on the mean (1SD) [cts]	1.66E+01	1.59E+01	2.42E+01	2.54E+01	3.42E+01	2.61E+01	5.31E+01	1.46E+02	2.89E+01	5.24E+01	3.85E+01	5.18E+01
Relative total combined uncertainty on the mean (1SD) [%]	0.62%	0.23%	0.28%	0.30%	0.38%	0.10%	0.16%	0.51%	0.32%	0.20%	0.12%	0.18%

- 0.01 s dwell time; 15  $\mu\text{m}$  spot size; 3.5  $\text{Jcm}^{-2}$
- 10 lines of 30s each @ different lasing parameters:
  - Test 1(light green) – 10 Hz, total gas flow of 0.65 LPM (0.35 LPM in MFC1 and 0.3 LPM in MFC2)
  - Test 2 (light blue) – 50 Hz, total flow of 0.75 LPM (0.45 LPM in MFC1 and 0.3 LPM in MFC2)
  - Test 3 (light orange) – 50 Hz, total flow of 1.05 LPM (0.525 LPM in MFC1 and 0.525 LPM in MFC2)

# Summary

- Laser ablation provides outstanding sample flexibility
- Capable of single spot or bulk analysis
- Capable of generating sample maps in either 2D or 3D
- Technique developed need for automated data interrogation software
- An excellent addition to your “analytical toolbelt”



## ADVANCED ANALYTICAL SCANNING TRANSMISSION ELECTRON MICROSCOPY (STEM), AND FUTURE DIRECTIONS

---

SPEAKER:

**Patrick Phillips, PhD**

Asst. TEM Product Manager,  
JEOL USA

March 24, 2022 | 11am PT



**COVALENT  
ACADEMY**

Advancements in  
Instrumentation Series

**Episode 31**

Want to learn more about Covalent's  
**Laser Ablation ICP-MS Services?**

**Talk with a Covalent Expert!**

Schedule your Appointment Now with Calendly

- link is in the chat -

## Science Forward.

Covalent delivers quality data and expert analysis for advanced materials and device innovation.


[All Analytical Services](#) 

 [Learn more about Advanced Search Beta](#)

[Contact Us](#)

## Access Covalent Portals

### Customer Access to Data & Community Content



The DATA PORTAL is used by Customers and Lab Partners for uploading and downloading data. It requires two-factor authentication and advanced password protection. Data Portal users have complete access through their home page on the portal to all Community content, and do not require a separate Community account.

LOGIN TO DATA PORTAL

### Community Content Access for **All Other Users**



The COVALENT COMMUNITY PORTAL requires password entry. It contains webinar and other metrology and characterization-related content that we believe would be useful and educational for the materials science innovation community. It does not provide access to any customer data and should only be used by individuals that are not Covalent customers or lab partners.

LOGIN TO COMMUNITY PORTAL



# Q & A Session



**COVALENT**  
**METROLOGY**

Thank you