

BATTERY MATERIAL QUALITY: THE IMPORTANCE OF PHYSICAL PROPERTIES

Anthony Chalou

Head of Market Development,
Anton Paar

June 20, 2024 at 8 AM Pacific Time

A background image showing a close-up of several red and silver cylindrical batteries. One battery in the foreground is in sharp focus, showing a white plus sign on its side. The background is blurred, showing more batteries and a bokeh effect of light.

**COVALENT
ACADEMY**

Industrial Applications of
Advanced Metrology
Episode 37



Modern, digitally-empowered analytical services platform delivering quality data and expert analysis to accelerate advanced materials and device innovation.



Comprehensive Solutions Stack

50+ Cutting-edge instruments in-house, 150+ Techniques

Analytical Services

Advanced Modeling

Method Development

Custom Consulting Solutions



High-touch, High-Quality Services

Ionic Membership Program

Enterprise Metrology Solutions

Instant, Secure Access to Data and Reports

Expanding Toolkit in Custom Digital Platform



Flexible Business Models

LiveView™ (real-time collaboration)

Co-op and Tool-share Opportunities

Training and Certification on Instrumentation

Laboratory Audits




Rich Network of Partnerships

Partner to World's Leading Instrument Manufacturers and Labs

Expanding access to Advanced Instruments and Analysis Tools

Lab Connections and Applications Learning



Who We Are, Who We Serve

50+ People, 14 PhDs

Comprehensive, Modern Analytical Capabilities

\$20M + State-of-art Lab in Sunnyvale, CA

600+ Clients, 15-30 new clients / week

Covalent's Analytical Services & Technical Groups



Electron Microscopy

- S/TEM
with EDS; EELS; Electron Diffraction; SAED
- FIB-SEM & HR-SEM
with EDS; EBSD; 3D Tomography
- Lamella Preparation
incl. specialized lift-outs



Failure Analysis

- Root-Cause Failure Analysis
- DPA / Mechanical X-section
- Dye & Pry Test
- Hot Spot Detection
- Emission Microscopy
- NIR / IR Imaging
- EBIC / OBIC failure analysis



Microscopy & Profilometry

- Chromatic Aberration
- Digital Optical Microscopy
- Laser Scanning Confocal Microscopy
- White Light Interferometry
- Scanning Acoustic Microscopy (SAM)



Mechanical Testing

- AFM & Advanced AFM Modes (EFM, KPFM, MFM, PFM, PiFM)
- Nano-indent / Nano-scratch
- Rheometry / Viscosity
- DMA / TMA (bend/stretch/compression)
- Tensile testing



Analytical Chemistry

- Mass Spectroscopy: ICP-MS and LA-ICP-MS; GCMS
- ICP-OES / GDOES
- Raman
- NMR (solid / liquid + 1,2,3 nuclei)
- XPS, UPS, ISS
- SIMS, TOF-SIMS



Misc. Material Properties

- Thermal Analysis: DSC, TGA
- Surface Zeta Potential
- Porometry / Pycnometry
- Gas Adsorption / Chemisorption
- Foam Density / Skeletal Density / Tap Density
- Particle Analysis: DLS / ELS / size distribution / zeta potential



X-ray Characterization

- X-Ray Diffraction (XRD)
- X-Ray Reflectometry (XRR)
- Micro-computed X-ray Tomography (Micro-CT)
- 2D / 2.5D / 3D X-ray Inspection & X-ray Radiography
- ED-XRF / WD-XRF



Optical Characterization

- Fourier Transformed Infrared Spectroscopy (FTIR and ATR-FTIR)
- Spectral Ellipsometry & Advanced Optical Modeling
- UV-Vis-NIR Spectroscopy

Today's webinar is in partnership with



Anton Paar

- **Partnership with Anton Paar announced in May, 2020**
 - Established the **Anton Paar Demonstration Facility in Covalent's Silicon Valley Laboratory** with goals of
 - Expanding industry access
 - Developing new analytical applications
 - Later expanded partnership to deliver industry-leading porous materials and powders analysis
 - Partners **continue to collaborate in advanced applications development**
- **Anton Paar Instruments** at Covalent Metrology include:
 - **SurPASS 3**
 - **Litesizer 500**
 - MCR 702 Rheometer / DMA
 - STeP 6 Nanoindentation platform
 - Ultrapyc 5000 Micro
 - Autosorb iQ C-XR-XR with CryoSync accessory
 - Porometer 3G and DualAutoTap
 - NEW Upgraded Nova 800 BET (Gas Adsorption) Analyzer

Other Covalent Partners



Rigaku



HORIBA
Scientific

ThermoFisher
SCIENTIFIC



KEYENCE



PVA TePla



Toray Research Center, Inc.

Anthony Chalou

*Head of Market Development,
Anton Paar*



- Anthony Chalou is the **Head of Market Development** at Anton Paar
- He holds **two Masters degrees in Applied Analytical Chemistry** and has **more than 20 years of experience with analytical instrumentation**, both in scientific applications as well as in sales and marketing
- His **current focus is on the characterization of physical properties of materials used in the R&D and QC of lithium-ion batteries**





Anton Paar


ANTON PAAR SOLUTIONS FOR
BATTERY MATERIAL CHARACTERIZATION

Anthony Chalou
Head of Market Development
Anton Paar GmbH

OUTLINE

- › Main components
- › LIB Types
- › LIB Chemistries
- › Cell formats
- › Pros & cons of common formats
- › Battery trilemma
- › Emerging trends
- › Physical parameters of raw materials
- › Physical parameters of slurries
- › Questions

MAIN COMPONENTS OF A LITHIUM ION BATTERY



TYPES OF LITHIUM ION BATTERIES



Ternary lithium ion (LIB)

- Most researched & used type of Li ion batteries
- Nickel-rich cathodes provide higher energy densities, especially when combined with Silicon based anodes
- Need reduction in production costs



Lithium Sulfur (Li-S)

- Sulfur is the main low-cost material used
- During discharge, lithium sulphides are formed at the cathode
- Low cycle life and safety concerns need to be further evaluated prior to full adoption

Estimated energy density range (Wh/kg)



Lithium Iron Phosphate (LFP)

- Mainly produced in China; few western manufacturers
- Cost advantage: No Co or Ni, and low cost of Iron and Phosphate
- Excellent cycle life, and improved safety
- Suffers from low energy density



Solid state lithium (SSB)

- The electrolyte is in solid form
- Very high energy density as well as safety, which make this type adoptable by the automotive industry
- Challenges in material compatibility and manufacturability in series are still issues that stand in the way of full market maturity

LITHIUM-BASED RECHARGEABLE SYSTEMS


| Technology | Anode | Cathode |
|---|-----------------|--------------------------|
| Lithium Cobalt Oxide 3.6V | Graphite –Gr/Si | LCO |
| Lithium Manganese Oxide Spinel 3.8V | Graphite –Gr/Si | LMO |
| Lithium Nickel Manganese Cobalt 3.7-3.8V | Graphite –Gr/Si | LNMC (111,532, 622, 811) |
| Lithium Nickel Cobalt Aluminum 3.65V | Graphite –Gr/Si | LNCA |
| Lithium Iron Phosphate 3.2V | Graphite –Gr/Si | LFP* |
| Lithium Manganese Iron Phosphate 3.7V | Graphite –Gr/Si | LMFP* |
| Lithium Titanite Oxide 2.3/2.4V | LTO | NMC, LMO, LFP |
| Lithium Nickel-Rich Oxide 3.8V | Graphite –Gr/Si | eLNO* |
| Lithium Metal 3.7V | Lithium Metal | NMC |
| <i>Lithium Nickel Manganese Oxide Spinel 3.8V (Mn Rich)**</i> | Graphite –Gr/Si | LNMO* |
| <i>Lithium Sulfur 2.4V**</i> | Lithium Metal | Sulfur* |

* Cobalt free

** R&D stage


MOST COMMON FORMATS OF BATTERY CELLS

Cylindrical



A spirally wound design (jelly-roll). Designated by size, e.g. 18650 cylindrical battery (Diameter: 18.6 mm, length: 65.2 mm; code for cylindrical shape: 0)

Prismatic




A prismatic design indicate a flat battery design. The stacks can be wound (as shown in the photo) or stacked (with alternating cathode/separator/anode structure). The stacks are usually inserted into rigid casing to form prismatic containers

Pouch



Rather than rigid metallic casing, conductive foil-tabs are welded to the electrodes and seal the battery fully. The stacks inside can be wound or stacked. Swelling and gassing could be a concern for pouch cells.

Coin cell




Also known as button cell. The cells are stacked into a tube. Most coin cells are single-use.

Source: J.M. Tarascon, Nature 2001, Sanyo and Panasonic

PROS AND CONS OF THE DIFFERENT FORMATS OF BATTERY CELLS

| Cylindrical | | Prismatic | | Pouch | |
|-----------------------------------|---|--|--------------------------------------|-----------------------------------|-----------------------------------|
| + | - | + | - | + | - |
| High energy density | Space inefficiency | High energy density | Mainly large and high-capacity cells | High volumetric energy density | Mechanical limitations |
| Wide operating voltage range | Low volumetric energy density (at pack level) | Wide operating voltage range | Requires special orders | Wide operating voltage range | Requires special orders |
| Wide power density range | Thick | Wide power density range | Thicker than pouch cells | Wide power density range | Swelling |
| Good cycle life (500-1000 cycles) | No custom sizes | Good cycle life (500-1000 cycles) | No standard sizing | Good cycle life (500-4000 cycles) | No internal safety devices |
| Low self-discharge (~1% / month) | Difficult thermal management | Low self-discharge (~1% / month) | Difficult thermal management | Low self-discharge (~1% / month) | Vulnerable to shock and vibration |
| Quick charge possible | | Quick charge possible | | Quick charge possible | |
| Internal safety devices | | Internal safety devices | | Better thermal management | |
| Standard sizes (no custom orders) | | Customization possible for mass production | | Easy to custom make | |
| Robust | | Robust | | Mass production for large cells | |
| Mass production | | Thin (Better pack volume) | | Thin (Better pack volume) | |

THE “IDEAL” BATTERY



EMERGING TRENDS



RECYCLING

- Black mass
- Hydrometallurgy
- Pyrometallurgy
- New regulations



SOLID STATE BATTERIES

- Enhanced safety
- Extended lifespan
- Increased energy density
- More expensive
- Challenges in scaling up



DRY COATING

- No solvent use
- Lower production cost
- Better performing batteries
- Challenges in production speed needs

BATTERY MATERIAL CHARACTERIZATION



Raw Material

- BET Surface area
- Pore size
- Pore volume
- Particle size
- Particle shape
- Skeletal density
- Tapped density
- Powder Rheology
- Defect ratio
- Crystallinity



Slurries

- Flow behavior
- Slurry viscosity
- Zeta potential
- Density
- Percent solids concentration
- Solvent concentration
- Particle size
- Particle shape



Coating

- Flow behavior
- Slurry viscosity
- Sedimentation stability
- Mechanical properties



Drying

- Adhesion
- Mechanical properties



Calendering

- Hardness & elastic modulus
- Adhesion
- Mechanical properties



Stacking

- Separator pore size
- Separator pore volume
- Mechanical properties



Electrolytes

- Separator pore size
- Separator pore volume
- Electrolyte viscosity
- Electrolyte density
- Electrolyte flash point




Battery Cell


- In operando characterization
- In situ characterization




IMPORTANT PHYSICAL PARAMETERS OF BATTERY RAW MATERIALS



BET SURFACE AREA



Eleven-point BET plot for NMC




Eleven-point BET plot for graphite




NOVA Series

- › Surface area of electrode material directly influences charge transfer between active materials and electrolyte
- › Surface area is optimized for electrical performance and to minimize adverse chemical reactions
- › Specific surface area of electrode materials is determined by the cryogenic adsorption of a gas (such as N₂)
- › Method uses gas sorption data to quantify the number of molecules of adsorbate in a monolayer on the surface of the material

PORE SIZE AND PORE VOLUME



N₂ isotherm on 3 types of hierarchical carbon samples (mesopores)




CO₂ isotherm on 3 types of hierarchical carbon samples (micropores)



NOVA Series

Kaomi Software




Merged CO₂ (273 K) and N₂ (77 K) DFT pore size distributions (blue) and cumulative pore volumes (red) for three carbons

PARTICLE PARAMETERS*


| Size | | Size distribution | | Shape | |
|-----------------------------------|---------------------|------------------------------|-----------------------|---------------------------------|--------------------------|
| Small size | | Uniform | Broad | Single crystal structure | |
| Shorter Li-ion transport distance | Larger SSA | Better cycle performance | Higher energy density | Higher electrode packed density | Difficult to manufacture |
| Better rate performance | More side reactions | Higher production difficulty | Poor cell homogeneity | Better cycle performance | Higher cost |

*Adapted from J. Zhang et al. "Balancing particle properties for practical lithium-ion batteries", *Particuology*, Vol. 61, (2022), P. 18-29

PARTICLE SIZE



Intensity-weighted particle size distribution of carbon black



Particle size distribution of three different cathode materials

Litesizer

Dynamic light scattering (DLS)
Electrophoretic light scattering (ELS)



0.3 nm

40 nm

10 µm




2.5 mm

Particle Size Analyzer (PSA)

Laser diffraction (LD)

PARTICLE SHAPE

Measurement results



Particle shape and size distribution of NMC cathode material measured with a dynamic image analyzer using a liquid flow dispersion unit

Particle shape affects parameters such as:

- › Homogeneity of the slurry
- › Uniformity of the coating on current collectors


For electrode material, near-spherical shape is preferred



Litesizer DIA

SKELETAL DENSITY

- › Skeletal density, or true density, is measured by gas pycnometry
- › Greater skeletal density of electrode material → higher energy density of batteries



Increase of skeletal density with increasing thermal treatment of nickel oxide electrode material

Principle of gas pycnometry:

Instead of weighing the displaced liquid, gas volumes are assessed via pressure measurements before and after expansion into an empty volume




Ultracyc Series

TAPPED DENSITY


- › Tapped density, also known as bulk or packing density, is an important parameter in battery material characterization, and incorporates:
 - › Volume of the solid phase
 - › Volume of pores within particles
 - › Volume of the spaces between particles
- › For solid battery material, a greater tapped density produces more efficient energy management

Example: LiNiCoMnO_2 cathode material, mixed particle sizes results in greater tapped density and better electrochemical performance than with single particle size systems



Autotap

POWDER RHEOLOGY




Compressibility of Si and Si/C samples as a function of normal stress

- › Si/C is denser than Si (also confirmed by gas pycnometry)
- › Si almost completely lacks compressibility (densification)
- › Si powder readily densifies under its weight; it is very flowable => it acts as a flow-aid




**MCR 302e /
Powder shear cell**



Powder flow parameters of Si/C samples as a function of moisture

- › Small effect of moisture levels between 0 % - 50 %
- › For moisture levels >50 % reduction in flowability, due to formation of menisci between particles
- › Important information for storage and transportation conditions

DEFECT RATIO




Raman spectra of graphene powder (black) and a graphene film (red).

- › “Graphitic” or “G” band at 1580 cm⁻¹: sharp signal in ordered graphite compared to amorphous carbon
- › “Defect” or “D” band around 1350 cm⁻¹: signal associated with the breathing vibrations of the carbon atoms, becomes visible around defects
- › The overtone of the D band, commonly referred to as the “2D” band, around 2680 cm⁻¹: for symmetry reasons, the 2D band in pure single-layer graphene is a perfect single-component peak
- › “Defect ratio” is the ratio of the D and G band intensities (I_D / I_G) that changes with the introduction of defects. Another ratio is the one between the G and the 2D band intensities (I_{2D} / I_G) which changes with the number of layers



Cora 5001

DEGREE OF CRYSTALLINITY



Diffractogram of a graphite anode material

- › Powders with high degree of crystallinity facilitate electron mobility in the battery cell
- › XRD gives information on the crystal structure, crystallite size, and phase purity of battery raw materials
- › XRD is capable of detecting minor impurities when comparing the diffractogram of a known material with that of a contaminated one
- › Specialized gas-tight sample holders also make it possible to measure air-sensitive materials without sample degradation

| Crystallite size (nm) | Degree of graphitization (%) |
|-----------------------|------------------------------|
| 41.2 | 92.7 |




Sample holders for air-sensitive battery materials



XRDynamic 500


IMPORTANT PHYSICAL PARAMETERS OF BATTERY SLURRIES



FLOW BEHAVIOR

| Speed [rpm] | Torque [%] | Temperature [°C] | Viscosity [mPa·s] |
|-------------|------------|------------------|-------------------|
| 60 | 82.4 | 40 | 5945 |

Dynamic viscosity of slurry with single point measurement





Calculated Herschel-Bulkley parameters to determine the flow behavior of slurry

The flow index is a description of the Newtonian or non-Newtonian behavior of the sample. A flow index that is smaller than 1 indicates that this slurry is shear-thinning – as is the case here for this slurry.




ViscoQC 300

SLURRY VISCOSITY



Viscosity curves measured after stepwise addition of components, representing the processing steps for mixing procedures A and B




Comparison of mixing procedures A and B in terms of total solid fraction and viscosity, at constant shear rate



MCR 302e

ZETA POTENTIAL (ζ)

- › The zeta potential describes the stability of colloidal dispersions.
- › Zeta potential is an important parameter for the characterization of particles in infusions and emulsions.
- › The higher the magnitude of the zeta potential (highly + OR highly -), the more stable the colloid.
- › Zeta potential is measured by electrophoretic light scattering (ELS)
- › During an ELS measurement an electric field is applied, which induces a collective motion of particles: particles with a higher zeta potential will move faster than less-charged particles.



**Omega
cuvette**



Litesizer 500
Dynamic light scattering (DLS)
Electrophoretic light scattering (ELS)

DENSITY AND PERCENT SOLIDS

- › Gas pycnometry can be used to determine the density and percent solids of electrode slurries
- › This is relevant for optimizing the surface area of the solid electrode




Ultrapyc Series

| % solids formulated | Measured density (g/cm ³) | % solids calculated |
|---------------------|---------------------------------------|---------------------|
| 41.45 | 1.3459 | 41.52 |
| 33.9 | 1.2682 | 34.27 |
| 26.30 | 1.1977 | 26.86 |
| 17.75 | 1.1226 | 17.96 |
| 15.22 | 1.1001 | 15.05 |
| 10.57 | 1.0688 | 10.81 |
| 6.41 | 1.0388 | 6.50 |
| 2.52 | 1.0121 | 2.45 |

Formulated and calculated % solids in slurries with measured densities

SOLVENT CONCENTRATION

- › N-Methyl-2-pyrrolidone (NMP) is a strongly polar aprotic solvent commonly used in battery slurry preparation to dissolve binders like PVDF
- › NMP is also used as a solvent for lithium salts in electrolytes
- › Because of its relative high price as well as its negative environmental impact, NMP is often recovered after use in the manufacturing process of lithium batteries
- › Refractive index correlates with the concentration of NMP, so it is a quick and inexpensive way to measure it from 0.0% to 99.9%




Refractive Index of different NMP concentrations measured at 20 °C and 589 nm




Abbemat 550

SEDIMENTATION STABILITY



Frequency sweeps of sample battery slurry A and battery slurry B.



- › A slurry is considered stable when the storage modulus G' is higher than the loss modulus G''
- › High frequencies enable evaluation of the samples' short-term behavior, while low frequencies simulate the behavior of the sample over a long time
- › $G'' > G'$ would be considered the cross-over point where the slurry sediments
- › In this example, slurry B exhibits better stability (cross-over at ~ 1 rad/s) compared to slurry A (cross-over at ~ 10 rad/s)



MCR 302e

QUESTIONS?



www.anton-paar.com/us-en/products/industries/application/lithium-ion-batteries/

View On-Demand Recordings in the Covalent Academy

Search resource topics

SEARCH

From
Select

To
Select

Type
Select

Duration
Select

Learning Center

- [Browse All Topics](#)

Register and start exploring at:
academy.covalentmetrology.com

PHOTO-INDUCED
FORCE MICROSCOPY
(PIFM): AUGMENTING
SURFACE ANALYSIS
WITH AFM CHEMICAL
MAPPING



Webinar

50 min

SURFACE CHARGE ON
COLLOIDS, AND
BEYOND: THE
COMPLEMENTARITY
OF SOLID- AND
SOLUTION-STATE
ZETA POTENTIAL
MEASUREMENT

Thomas Luxbacher, PhD
Principal Scientist,
Anton Paar

March 9, 2023 | 11am PT



COVALENT
ACADEMY
Industrial Applications of
Advanced Metrology
Episode 34

Webinar

60 min

NANOSIMS:
HIGH SENSITIVITY
IMAGING ANALYSIS
FOR DEVICES

SPEAKER:
Junichiro Sameshima, PhD
Senior Manager of
Surface Science Laboratory,
Toray Research Center, Inc.
September 22, 2022 | 11am PT

Toray Research Center, Inc.



COVALENT
ACADEMY
Advancements in
Instrumentation Series
Episode 33

Webinar

60 min

CHARACTERIZATION
OF CLIMATE
RESEARCH



ADVANCED ANALYTICAL
SCANNING
TRANSMISSION



LASER ABLATION
INDUCTIVELY
COUPLED PLASMA





Q & A Session



**COVALENT
METROLOGY**

Thank you.