



COVALENT
METROLOGY

Welcome

**ADVENTURES IN
WAVELENGTH
DISPERSIVE X-RAY
FLUORESCENCE
(WDXRF): FLEXIBLE
ELEMENT ANALYSIS
FOR THIN FILMS
AND MORE**

SPEAKER:

Meredith Beebe

Semiconductor X-ray
Metrology Specialist,
Rigaku

January 27, 2022 | 11am PT

 The Covalent Academy logo graphic is a circular emblem with a dark blue background and a yellow border. It contains the text 'COVALENT ACADEMY' in white, 'Advancements in Instrumentation Series' in a smaller white font, and 'Episode 29' in a bold white font. Four yellow dots are positioned at the top, bottom, left, and right of the inner circle.

**COVALENT
ACADEMY**


Advancements in
Instrumentation Series

Episode 29



COVALENT METROLOGY

Silicon Valley-based analytical labs and platform delivering quality data and expert analysis for advanced materials and device innovation



Comprehensive Solutions Stack

50+ cutting-edge instruments, offering 100+ Techniques

Analytical Services

Advanced Modeling

Method Development

Temp. Staffing Solutions



Affordable and Fast

Fast Turnaround Times, No Expedite Fees

Volume Savings

Instant Access to Data and Reports in Secure Portal



Flexible Business Model

Custom Consulting Solutions and Certified Onsite Support

Training and Certification on Instrumentation

Co-op and Tool-Share Opportunities

Laboratory Audits



Rich Network of Partnerships

Partner to World's Leading Instrument Manufacturers and Labs

Expanding Instrumentation, Lab Connections and Learning



Who We Are, Who We Serve

500+ Clients, 40-60 Added / Quarter

50+ People, 14 PhDs

Cutting-edge Analytical Capabilities

Lab Location: Sunnyvale, CA

Covalent Technical Groups and Organization

4

PCBA, Semiconductor, and Electronic Device Metrology & Failure Analysis

- DPA / Mechanical Cross-section
- Dye & Pry Test
- EBIC / OBIC failure analysis
- Hot Spot Detection
- IR Imaging / Emission Microscopy
- NIR Imaging
- Root-Cause Failure Analysis

Electron Microscopy and Scanning Probe Microscopy

- AFM & Advanced AFM Modes (EFM, KPFM, MFM, PFM)
- Scanning Acoustic Microscopy (SAM)
- SEM (+ EDS)
- FIB-SEM (+ EDS)
- S/TEM (+ EDS / + EELS)
- Nano-indent / Nano-scratch

Optical Microscopy & Spectroscopy

- Chromatic Aberration
- Digital Optical Microscopy
- FTIR and ATR-FTIR
- Laser Scanning Confocal Microscopy
- Spectral Ellipsometry
- UV-Vis-NIR Spectroscopy
- White Light Interferometry

X-Ray Characterization

- X-Ray Diffraction (XRD)
- X-Ray Reflectometry (XRR)
- Micron-spot ED-XRF
- WDXRF
- Micro-computed X-ray Tomography (Micro-CT)
- 2D X-ray Inspection & X-ray Radiography

Elemental / Chemical Composition Analysis

- EPMA
- GD-OES
- GC-MS
- ICP-MS and LA-ICP-MS
- Raman Microscopy & Spectroscopy
- NMR (1D or 2D; solid / liquid)

Particle Analysis

- Dynamic Light Scattering (DLS)
- Laser Diffraction Particle Size Analysis (PSA)
- Particle Zeta Potential

Material Property Characterization

- DSC
- DMA & TMA
- Rheometry
- TGA
- Surface Zeta Potential
- Porometry / Porosity
- Gas Adsorption
- Gas Pycnometry
- Foam Density
- Tap Density

Surface Spectroscopy Analysis

- Dynamic-SIMS
- ToF-SIMS (Static-SIMS)
- Ion Scattering Spectroscopy (ISS)
- Ultraviolet Photoelectron Spectroscopy (UPS)
- X-ray Photoelectron Spectroscopy (XPS)



- Partnership announced Fall, 2018
- Covalent is home to Rigaku's **Silicon Valley Semiconductor Division Metrology Lab:**
 - SmartLab XRD / XRR with full set of optics + heated stage
 - CT Lab HX 130 X-ray Micro-CT
 - AZX 400 WDXRF with 200 mm and 300 mm cassettes / FOUP (in clean booth)

Other Covalent Partners



Meredith Beebe


Semiconductor X-ray Metrology Specialist,
Rigaku

- Over 20 years in semiconductor industry spanning roles both as an end user in the fab (as a process characterization engineer) and as an applications expert for multiple vendors
- Subject matter expert known across the industry for trace metal and contamination analyses employing X-rays, VPD-ICP-MS, SIMS, AFM, and Auger
- Provides a platform to address ever-changing metrology demands for yield improvement and advanced device development
- Helps Rigaku customers to understand and develop improved strategies for contamination monitoring, metrology optimization, and she provides key insights from direct customer experience to contribute to Rigaku's leadership in X-ray metrology




OUTLINE

- WDXRF technique overview
 - Principles of X-ray fluorescence
 - Sequential and simultaneous spectrometers
 - Qualitative and quantitative methods
 - Comparison of WDXRF to EDXRF
- Rigaku AZX 400 spectrometer overview
- Application studies
 - PZT
 - Sputtering target
 - MRAM
- Summary/Q&A




FLEXIBILITY OF RIGAKU'S AZX 400 WDXRF

- Non-destructive detection of a wide range of elements from Be-U
 - Example application– MRAM
 - Particularly suited to light elements such as boron through aluminum
 - Ultra-thin MgO film @ 1 nm
- Flexible sample handling supporting enterprise needs from process R&D to manufacturing
 - Example application– Sputtering target
 - Multiple sample adaptors allow for various shapes and sizes of coupons, targets, and wafer samples
 - Sample size range up to 380 mm
- Film thickness measurements from angstroms to microns
 - Example application– PZT
 - High-precision measurements



BASIC PRINCIPLES OF X-RAY FLUORESCENCE



X-ray (high-energy photon) is directed at the sample.

Incident X-ray with sufficient energy ejects an electron from an atom leaving a vacancy.




This vacancy puts the atom in an unstable excited state with a higher energy.

Atom relaxes to restore to its original configuration. It transfers an electron from an outer shell to the lower energy level to fill vacancy.

The replacing electron loses energy, which is radiated as X-rays of a characteristic energy.



The intensity of fluorescent X-ray photons is proportional to a given element.

FUNDAMENTALS OF WDXRF




Wavelength-energy relationship:

$$\lambda (\text{\AA}) = \frac{12.4}{E (\text{keV})}$$




SEQUENTIAL WDXRF




Sequential spectrometers use a moving detector on a goniometer to move through an angular range to measure the intensities of many different wavelengths.

SIMULTANEOUS WDXRF



Simultaneous spectrometers are equipped with a set of fixed detection systems, where each system measures the radiation of a specific element.


QUALITATIVE SCANS



Non-quantitative scans of elements across the periodic table can be set up to identify elements present in the sample

FILM THICKNESS AND X-RAY INTENSITY

Film thickness is obtained by calibration based on the relationship between thickness and X-ray intensity.




QUANTITATIVE METHODS


Standardization models to address thin film and bulk samples

- Fundamental parameter (FP)
 - Pure element or “type standards”– quantitative
 - Type standards are similar in composition to “unknowns” to be measured
 - Improves accuracy of FP calibration
 - Library standards– semi-quantitative
 - Computer-stored calibration parameters
 - Helpful when calibration standards are not available
 - Density is important for determining linear thickness
- Empirical models
 - Relies on a suite of calibration standards covering the range and thickness values expected of “unknown” product

FP AND EMPIRICAL CALIBRATIONS





- Theoretical vs measured intensity
- Theoretically linear
- Thickness and composition analysis can be done simultaneously with multiple layers
- Few standards are enough, pure metals can be used




- Measured intensity vs standard value
- Not usually linear
- Analyzes only thickness or composition
- Range of standards necessary for each element, thickness, and composition

WDXRF AND EDXRF COMPARISON

	WDXRF	EDXRF
System Set-up	 <p>The diagram shows an X-ray source emitting a beam through a slit onto an analyzing crystal. The crystal diffracts the beam, and a detector captures the resulting signal.</p>	 <p>The diagram shows an X-ray source emitting a beam directly onto a sample, with a detector positioned to capture the emitted fluorescence.</p>
Measurable Elements	Be - U	Al - U
Quantitative Method	<p>Advanced FP Method is available</p> <p>FP method enables thickness and composition analysis with one recipe.</p>	<p>Usually Empirical Method</p> <p>Empirical method cannot analyze thickness and composition simultaneously and requires a range of standards for thickness and composition calibration.</p>
Analyzable Structure	<p>Multiple stacked layers can be analyzed</p> <p>FP method considers absorption effects of other layers and can analyze complex compounds and multi-layered samples.</p>	<p>Single layer or A few layers of simple structure</p>
Analyzable Thickness	< 1 nm or more	10 nm or more

RESOLUTION COMPARISON OF WDXRF AND EDXRF


1 μm Al film



WDXRF can clearly resolve the Al-K α peak (1.487 keV) without influence from the substrate Si-K α peak (1.740 keV).

RESOLUTION COMPARISON OF WDXRF AND EDXRF

GST film with Ge, Sb, and Te



Clear separation of Sb-La (3.606 keV) and Te-La (3.770 keV) peaks in GST film with WDXRF.


SEQUENTIAL WAVELENGTH DISPERSIVE XRF RIGAKU AZX 400



- 4 kW sealed X-ray tube, sequential type goniometer, primary beam filter; Measurement spot sizes 30, 20, 10, 1, and 0.5 mm in diameter
- Flexibility to measure a variety of sample types, including 50 - 300 mm wafers, coupons, and sputtering targets. Auto wafer loader option available.
- Analyzes elements from Be to U with a range from ppm to % or thickness range of Å to μm
- Well-suited for process R&D and low-volume, high product mix environments
- SEMI and CE compliant with a small footprint, 1376 (W) x 1710 (H) x 890 (D) mm

3 DIFFERENT APPLICATIONS OF WDXRF


- MRAM
- Sputtering target
- PZT



APPLICATION 1: MRAM


Magnetoresistive random access memory

- MRAM is a non-volatile random-access memory which stores data in magnetic domains.
- Faster writing speeds and re-writing times than flash memory.
- MRAM basic structure is a magnetic tunnel junction (MTJ).



MTJ structure consists of two ferromagnetic (FM) layers separated by an insulating tunnel barrier.

APPLICATION 1: MRAM



Good adventure for WDXRF	
✓	MTJ structure has light elements of Mg and B.
✓	High elemental resolution necessary to separate overlapping peaks.
✓	High precision required for thickness and composition repeatability measurements


- Thickness and composition monitoring of each layer are necessary quality controls for manufacturing.

ANALYSIS OF MgO LAYER


Ta layer is often used as a cap and/or barrier around the MTJ structure.

In this case, the Ta-Ma peak that is near Mg-Ka (1.254 keV) is fully resolved.

With EDXRF, the Mg-Ka peak can overlap with the high background of Si-Ka (1.74 keV) and Ta-Ma (1.77 keV).



ANALYSIS RESULTS OF MgO LAYER




10 times measurement of each sample

	MgO Thickness			
Target Thickness	0.4	0.8	1.2	1.6
Unit	Nm			
Data No.	10	10	10	10
AVERAGE	0.40	0.80	1.20	1.60
MAX	0.40	0.81	1.21	1.61
MIN	0.39	0.80	1.20	1.60
RANGE	0.01	0.01	0.01	0.01
Std. Dev.	0.0035	0.0042	0.0046	0.0052
R.S.D. (%)	0.88	0.52	0.39	0.32

WDXRF FP calibration of ultra-thin Mg thickness range from 4-16 Å.
 This technique can see 1 Å steps and ultra-thin thicknesses of 5 Å, which can be difficult for other methods.

APPLICATION 2: SPUTTERING TARGET ANALYSIS


- Sputtering targets are used to create thin films in sputter deposition tools.
- Composition and trace element detection is important to manufacturers.
- Stoichiometry of target composition can change with use



Good adventure for WDXRF

✓	Non-destructive and no sample prep needed for analysis
✓	Wide elemental range (Be-U) for various sputtering target materials
✓	Library standards can be used if standards are not readily available
✓	AZX 400 can accommodate large target samples or wafer samples from sputter target deposition

TARGET QUALITATIVE SCAN



Huge peaks are observed for Pb, Zr, and Ti dominant elements.

AZX 400 primary filters are used in analysis to suppress the background and high X-ray intensity for these elements and avoid counting loss errors for Pb and Zr.

This target also had trace amounts of Si and Al.

TARGET REPEATABILITY MEASUREMENT

Component	PbO Comp.	ZrO2 Comp.	TiO2 Comp.	SiO2 Comp.	Al2O3 Comp.
Unit	mol%	mol%	mol%	mol%	mol%
Data No.	10	10	10	10	10
AVERAGE	54.798	23.407	21.659	0.087	0.050
MAX.	54.849	23.438	21.722	0.089	0.053
MIN.	54.758	23.373	21.575	0.082	0.048
RANGE	0.091	0.065	0.147	0.007	0.005
Std. Dev.	0.025	0.023	0.041	0.002	0.002
R.S.D.(%)	0.05	0.10	0.19	2.65	3.05

Table: 10 times measurement of center point


Element Spectrum	Pb La	Zr Ka	Ti Ka	Al Ka	Si Ka
kV-mA	50kV – 60mA			30kV – 100mA	
Primary Filter	*F-Cu	*F-Cu	*F-Sn	OUT	OUT
Diameter	20 mm				
Slit	S2	S2	S2	S4	S4
Analyzing Crystal	LiF1	LiF1	LiF1	PET	PET
Detector	SC	SC	SC	F-PC	F-PC
Counting Time (s) Peak	10	10	10	20	20
Background	5	5	5	10	10

*Usually, primary filters are used to reduce the diffraction interference, but for this measurement they were used to reduce background and high intensity.

25 PT MAPPING ANALYSIS

Component	PbO Comp.	ZrO2 Comp.	TiO2 Comp.	SiO2 Comp.	Al2O3 Comp.
Unit	mol%	mol%	mol%	mol%	mol%
Data No.	25	25	25	25	25
AVERAGE	54.782	23.427	21.654	0.085	0.053
MAX.	54.888	23.480	21.860	0.112	0.099
MIN.	54.620	23.353	21.501	0.065	0.040
RANGE	0.268	0.127	0.359	0.047	0.059
Std. Dev.	0.064	0.034	0.090	0.013	0.011
R.S.D.(%)	0.12	0.14	0.41	15.49	19.96




Table: 25-point mapping test



25 pt map

SiO2 (mol%)

Al2O3 (mol%)



Maps of the trace SiO2 (mol%) and Al2O3 (mol%) detected

APPLICATION 3: PZT, LEAD ZIRCONATE TITANATE



- Piezo-electric material is in various sensors and MEMS devices in ink-jet printers, gyroscopes, mirrors for wearable devices, and more.
- Both the PZT thickness and composition are critical to device performance.

PZT (PbO , ZrO_2 , TiO_2)
Pt
Ti (TiO_x): Fixed Value
SiO_2 : Fixed Value
Si substrate

Good adventure for WDXRF	
✓	High spectral resolution needed to separate Pb, Pt, Ti overlapping peaks
✓	FP analysis can analyze complicated multilayer stacks for film thickness and composition
✓	AZX 400 can accommodate PZT wafer samples, coupons, and pure metal disc standards

PZT ANALYSIS WITH WDXRF


- FP analysis allows for TiOx layer thickness to be a fixed value.

This enables Ti to be measured and accounted for from the PZT layer.

- Pure metal standards of Pb, Zr, and Ti can also be used.

PZT (PbO, ZrO ₂ , TiO ₂)
Pt
Ti (TiO _x): Fixed Value
SiO ₂ : Fixed Value
Si substrate

QUALITATIVE CHART OF PZT SAMPLE



1 μ m PZT and 100 nm Pt spectra.
WDXRF high spectral resolution resolves the Pb-L line, Pt-L line, and Zr-K line that overlap each other in EDXRF.

PZT ANALYSIS RESULTS

Layer	PZT						Pt
Component	Thk.	PbO Conc.	ZrO2 Conc.	TiO2 Conc.	Pb / (Zr+Ti)	Zr / (Zr+Ti)	Thk.
Unit	nm	mol%	mol%	mol%	-	-	nm
AVERAGE	3282.1	49.42	26.27	24.31	0.977	0.519	100.0
MAX.	3288.7	49.48	26.35	24.38	0.979	0.521	100.2
MIN.	3272.3	49.37	26.19	24.23	0.975	0.518	99.7
RANGE	16.4	0.11	0.15	0.15	0.004	0.003	0.5
Std. Dev.	4.7	0.034	0.041	0.052	0.0013	0.0009	0.15
R.S.D.(%)	0.14	0.07	0.16	0.21	0.14	0.17	0.15

Table: 10 times repeatability measurement of PZT/Pt multilayer

*Measurement from fixed channel WDXRF system

Thickness and composition of the PZT and Pt layer are calculated.



10x repeat thickness measurements <.2% RSDs.

Composition of PbO, ZrO₂, TiO₂ are calculated.

Pb/(Zr+Ti) and Zr/(Zr+Ti) composition ratios are also determined from measured mol%.

AZX 400 MEASUREMENT FLEXIBILITY

- Different sample adaptors to accommodate a range of sample sizes
- Vacuum option
- Library standards are available for semi-quantitative analysis if standards are not readily available.
- FP method and high peak resolution enable resolution of complicated spectra and overlapping peaks.
- Primary beam filters – **diffraction elimination**



WDXRF THIN FILM APPLICATIONS

- **Semiconductor films**

Doped poly-Si (B, P, C, N, As), SiGe

- **Insulative films**

Insulator film:

PSG, BPSG

Low-k film:

AsSG, FSG SiOC

Nitride film:

Si₃N₄, SiON

High-k film:

La₂O₃, HfO₂, Ta₂O₅, Al₂O₃

- **Metal, conductive films**

Metal alloy film:

AlSiCu, AlCu, TiW, TaAl

Metal film:

W, Mo, Ti, Co, Al, Cu, Ir, Pt, Ru

Silicide film:

WSix, MoSix, TiSix, CoSix, NiSix

Nitride film (Barrier):

TiN, TaN, WN

- **Other**

Power Device:

Au(Ag)/Ni/Ti/Al (backside electrode)

Universal Memory:



MRAM (MgO, CoFeB), GST, PZT

MEMS (piezo, mobile):

PZT, (Sc,Ta)AlN, SAW, F-BAR


Other Purposes:

LED, sputtering target, photomask, F, Cl residue, C film



KEY FEATURES OF WDXRF

- Non-destructive detection for wide range elements from
- Ultra thin film thickness measurements from 0.1 nm
- High peak resolution and fundamental parameter (FP) analysis can analyze multilayer stacks with complicated compounds
- Multiple sample adaptors allow for various shapes, coupons, and wafer samples from 50 mm- 380 mm
- Library standards in software provides semi-quantitative measurement analysis without standards
- Qualitative measurements scan for smaller element groups or from Be-U



LASER ABLATION INDUCTIVELY COUPLED PLASMA MASS SPECTROSCOPY: NOT JUST ROCKS

SPEAKER:

Lucas Smith

Director of Business Development
for the Americas,
Teledyne CETAC

February 10, 2022 | 11am PT



COVALENT ACADEMY

Advancements in
Instrumentation Series

Episode 30



Want to learn more about Covalent's
WDXRF and Spectroscopy Services?

Talk with a Covalent Expert!

Schedule your Appointment Now with Calendly

- link is in the chat -

Science Forward.

Covalent delivers quality data and expert analysis for advanced materials and device innovation.


[All Analytical Services](#) 

 [Learn more about Advanced Search Beta](#)

[Contact Us](#)

Access Covalent Portals

Customer Access to Data & Community Content



The DATA PORTAL is used by Customers and Lab Partners for uploading and downloading data. It requires two-factor authentication and advanced password protection. Data Portal users have complete access through their home page on the portal to all Community content, and do not require a separate Community account.

[LOGIN TO DATA PORTAL](#)

Community Content Access for **All Other Users**



The COVALENT COMMUNITY PORTAL requires password entry. It contains webinar and other metrology and characterization-related content that we believe would be useful and educational for the materials science innovation community. It does not provide access to any customer data and should only be used by individuals that are not Covalent customers or lab partners.

[LOGIN TO COMMUNITY PORTAL](#)



Q & A Session



COVALENT
METROLOGY

Thank you