



**COVALENT
METROLOGY**

Welcome

CHARACTERIZATION OF CLIMATE BENEFICIAL MATERIALS BY GAS SORPTION

SPEAKER:

Martin Thomas, PhD

Lead Scientist,
Anton Paar Quantatec

May 19, 2022 | 11am PT




COVALENT ACADEMY

Advancements in
Instrumentation Series

Episode 32



Silicon Valley-based analytical labs and platform delivering quality data and expert analysis for advanced materials and device innovation



Comprehensive Solutions Stack

50+ cutting-edge instruments, offering 100+ techniques

Analytical Services

Advanced Modeling

Method Development

Temporary Staffing Solutions



Affordable and Fast

Fast Turnaround Times, No Expedite Fees

Volume Savings

Instant Access to Data and Reports in Secure Portal



Flexible Business Model

Custom Consulting Solutions and Certified On-site Support

Training and Certification on Instrumentation

Co-op and Tool-Share Opportunities


Laboratory Audits



Rich Network of Partnerships

Partner to the World's Leading Instrument Manufacturers and Labs

Expanding Instrumentation, Lab Connections, and Learning



Who We Are, Who We Serve

500+ Clients, 40 – 60 Added / Quarter

60+ Employees, 14 PhDs

Cutting-edge Analytical Capabilities

Silicon Valley Lab located in Sunnyvale, CA

Covalent Technical Groups & Organization



PCBA, Semiconductor, and Electronic Device Metrology & Failure Analysis

- DPA / Mechanical Cross-section
- Dye & Pry Test
- EBIC / OBIC failure analysis
- Hot Spot Detection
- IR Imaging / Emission Microscopy
- NIR Imaging
- Root-Cause Failure Analysis

Electron Microscopy and Scanning Probe Microscopy

- AFM & Advanced AFM Modes (EFM, KPFM, MFM, PFM)
- Scanning Acoustic Microscopy (SAM)
- SEM (+ EDS)
- FIB-SEM (+ EDS)
- S/TEM (+ EDS / + EELS)
- Nano-indent / Nano-scratch

Optical Microscopy & Spectroscopy

- Chromatic Aberration
- Digital Optical Microscopy
- FTIR and ATR-FTIR
- Laser Scanning Confocal Microscopy
- Spectral Ellipsometry
- UV-Vis-NIR Spectroscopy
- White Light Interferometry

X-Ray Characterization

- X-Ray Diffraction (XRD)
- X-Ray Reflectometry (XRR)
- Micron-spot ED-XRF
- WDXRF
- Micro-computed X-ray Tomography (Micro-CT)
- 2D X-ray Inspection & X-ray Radiography

Elemental / Chemical Composition Analysis

- EPMA
- GD-OES
- GC-MS
- ICP-MS and LA-ICP-MS
- Raman Microscopy & Spectroscopy
- NMR (1D or 2D; solid / liquid)

Particle Analysis

- Dynamic Light Scattering (DLS)
- Laser Diffraction Particle Size Analysis (PSA)
- Particle Zeta Potential

Material Property Characterization

- DSC
- DMA & TMA
- Rheometry
- TGA
- Surface Zeta Potential
- Porometry / Porosity
- Gas Adsorption
- Gas Pycnometry
- Foam Density
- Tap Density

Surface Spectroscopy Analysis

- Dynamic-SIMS
- ToF-SIMS (Static-SIMS)
- Ion Scattering Spectroscopy (ISS)
- Ultraviolet Photoelectron Spectroscopy (UPS)
- X-ray Photoelectron Spectroscopy (XPS)



- **Partnership with Anton Paar initially announced in May, 2020**
 - Established new Anton Paar Demonstration Facility in Covalent's Silicon Valley Laboratory to expand industry access and develop new analytical applications
 - **Initial instruments installed include:** MCR702 Rheometer / DMA, SurPASS 3, and Litesizer 500; STeP 6 Nanoindentation platform added in Fall 2021
- Anton Paar and Covalent later expanded partnership to **deliver industry-leading porous materials and powders analysis throughout the West Coast**
 - Anton Paar acquired Quantachrome in 2018 and in 2021 added new instrumentation to the shared Demo Facility at Covalent
 - **Instruments include:** Autosorb iQ C-XR-XR with CryoSync accessory, Ultrapyc 5000 Micro, Porometer 3G, and DualAutoTap system
 - **NEW Upgraded Nova 800 BET (Gas Adsorption) Analyzer** - for true surface area measurement and pore size distribution

Other Covalent Partners



Martin Thomas, PhD

Lead Scientist

Anton Paar Quantachrome

- Product Manager
- 30 Years with Quantachrome QuantaTec
- Obtained 10 patents in the field of porous materials analysis
- Co-author of “*Characterization of Porous Solids and Powders Surface Area, Pore Size, and Density*” (Springer)



CHARACTERIZATION OF CLIMATE BENEFICIAL MATERIALS BY GAS SORPTION

WHICH MATERIAL PROPERTIES CAN
GAS SORPTION BE USED TO
MEASURE?

WHY IS **GAS SORPTION** ANALYSIS
ADVANTAGEOUS FOR PRODUCT
RESEARCH RELEVANT TO CLIMATE
CHANGE AND SUSTAINABILITY?


HOW DOES **GAS SORPTION**
ANALYSIS WORK?

WHAT ARE SOME ADVANCED MODES OF
GAS SORPTION ANALYSIS THAT HAVE
EMERGED RECENTLY?




WHICH MATERIAL PROPERTIES CAN GAS SORPTION BE USED TO MEASURE?

- › Surface Area
- › Pore Volume
- › Pore Size Distribution
- › Absorption Capacity



FOR WHAT TYPES OF MATERIALS?

- › **Battery** cathode and anode powders
- › **Adsorbents** for air and water purification
- › **Catalysts** for pollution control
- › **Greenhouse gas storage** materials
- › **Hydrogen storage** materials



WHY IS GAS SORPTION ANALYSIS ADVANTAGEOUS FOR PRODUCT RESEARCH RELEVANT TO CLIMATE CHANGE AND SUSTAINABILITY?

- › **Surface area** is the quantitative expression of the extent of the interface between a solid and its surrounds... be it an adsorber of a toxic gas, or the electrodes of a lithium battery

Gas purification/separation, battery performance



WHY IS GAS SORPTION ANALYSIS ADVANTAGEOUS FOR PRODUCT RESEARCH RELEVANT TO CLIMATE CHANGE AND SUSTAINABILITY?

- › Surface area can further be qualified as **catalytically active surface area** – essentially exactly how efficient is your emission control catalyst.

Mitigating automotive/industrial emissions



WHY IS GAS SORPTION ANALYSIS ADVANTAGEOUS FOR PRODUCT RESEARCH RELEVANT TO CLIMATE CHANGE AND SUSTAINABILITY?

- › **Pore volume** represents the adsorption capacity primarily in high pressure applications where adsorption is not just 2-dimensional (limited to the surface) but 3-dimensional!

Greenhouse gas sequestration and storage



CARBON CAPTURE AND STORAGE


Functional solids for adsorption of CO_2 and CH_4 :

Geostores:

Oil/gas bearing rock (shale)
Coal

Synthetic carbons

Metal organic frameworks (MOFs)



WHY IS GAS SORPTION ANALYSIS ADVANTAGEOUS FOR PRODUCT RESEARCH RELEVANT TO CLIMATE CHANGE AND SUSTAINABILITY?

- › **Pore size** distribution controls diffusion of fluids and reactants in and out of the pore structure, and can impart a molecular sieving effect.

Minimizing rate limiting diffusion



WHY IS GAS SORPTION ANALYSIS ADVANTAGEOUS FOR PRODUCT RESEARCH RELEVANT TO CLIMATE CHANGE AND SUSTAINABILITY?

- > Absorption tells us how much hydrogen, for example, can be reversibly stored as a chemical compound in certain hydride forming metals.

On demand hydrogen generation for mobile applications




**CRYOGENIC ANALYSIS FOR
CLEAN **BATTERY** MATERIALS**

**REACTIVE GAS CHEMISORPTION TO
ANALYZE **ENVIRONMENTAL CATALYSTS****

**HIGH PRESSURE GAS SORPTION TO
PROBE **GREENHOUSE GAS**
SEQUESTRATION MATERIALS...**

...AND **FUEL CELL AND HYDROGEN
STORAGE MATERIALS**






GAS SORPTION ESSENTIALS


HOW DOES GAS SORPTION WORK?

- > **Physisorption**
(inert gases, low temperatures)
- > **Chemisorption**
(reactive gases, moderate temperatures)
- > **Absorption**
(solubility, bulk compound formation)
- > **Vacuum-volumetric** apparatus
- > Pressure
- > Temperature




WHAT IS PHYSISORPTION?

- > A weak & reversible physical interaction between a solid and an adsorbate (gas or vapor).
- > Used to measure surface area, pore size and pore volume.




WHAT IS CHEMISORPTION?

- › A strong and largely irreversible reaction between a solid and an adsorbate (gas or vapor).
- › Used to measure the extent of catalytically active area




WHAT IS HIGH PRESSURE SORPTION?

- › Can be physisorption, can be chemical absorption (bulk compound)
- › High pressures needed to overcome super-ambient or supercritical temperatures
- › Usually corresponds to process/storage temperatures



RECOMMENDATIONS

Material	Application	Property	Method	Gas
Mixed metal oxides	Battery cathode	Surface area	Physisorption	N ₂ @ 77K
Graphite	Battery anode	Surface area	Physisorption	N ₂ @ 77K
Activated carbon, zeolites	Gas/water purification	Surface area and micropore size distribution	Physisorption	N ₂ @ 77K Ar @ 87K CO ₂ @ 273K
Supported metal catalysts	Emission control, pollution control	Active metal area	Chemisorption	H ₂ , CO ~300K
Coal, shale	Carbon(CO ₂ /CH ₄) sequestration	Micropore volume	High pressure sorption	CO ₂ /CH ₄ >300K
Activated carbons	Hydrogen storage	Micropore volume	High pressure sorption	H ₂ ≤ 300K
Metals/alloys	Hydrogen storage	Hydride formation	High pressure sorption	H ₂ ~300K




SURFACE AREA AND PORE SIZE ESSENTIALS



SURFACE AREA

- Surface Area is the **extent of the interface** through which a solid quantitatively interacts
- Surface area is **created** by division of particles (size reduction) and the generation of porosity (drying, decomposition).

“low area”




SURFACE AREA BATTERIES, BATTERY MINERALS

- › More area = more interfacial current
- › More area = more side reactions and ageing

$$S = \frac{6}{\rho D} ?$$

Surface area S of spherical particles calculated from diameter D




*Realistic surface area from
gas adsorption*



INTERNAL/EXTERNAL SURFACE AREA

Batteries, Battery minerals


macroporous $\text{LiNi}_{0.5}\text{Mn}_{0.5}\text{O}_2$ ($8.87 \text{ m}^2\text{g}^{-1}$) from solution vs solid state reaction method ($2.93 \text{ m}^2\text{g}^{-1}$)



CATHODE SURFACE AREA




- › Higher surface area significantly reduces thermal stability
- › Activation energy of cathode decomposition decreases with increasing surface area
- › High surface area materials react more exothermally with electrolyte



EXPERIMENTAL DATA – CALCULATION

MODEL – B.E.T. AREA


> Isotherm > monolayer > area



Brunauer, Stephen; Emmett, P. H.; Teller, Edward (1938). *J. Am. Chem. Soc.* **60** (2): 309–319.

WHY PORE SIZE DISTRIBUTIONS?


= Pore volume and how it is distributed amongst pores of different size



- Pores **create** surface area
- Pore size and network **controls** fluid migration e.g. *reactants in and products out of a catalyst particle*
- Smaller pores **decrease** diffusion kinetics
- Smallest pores restrict access completely

PORE SIZE AND VOLUME

Are expressed both as cumulative and derivative pore size distributions





PORE SIZE

Catalyst supports

- › Silica (MS-1)
- › Alumina (MS-2)
- › Titania (MS-3)

Nitrogen, 77.4K	MS-1	MS-2	MS-3
BET surface area (m ² /g)	662	211	98
Total pore volume (mL/g)	0.705	0.376	0.235

Mesoporous Catalyst Supports – Nova x00 Series,
Anton Paar Application Report (2022)



ACTIVATED CARBON

- > Micropores
- > CO₂, 273K
- > NLDFT
- > Pore size merge




Figure 1: CO₂ (273 K) isotherm on activated carbon







Figure 2: NLDFT pore size distribution calculated from the CO₂ (273 K) isotherm


Activated Carbon – Micropore Size Distribution using CO₂ Adsorption – Nova x00 Series, Anton Paar Material Brief (2022)

SUPERCAPACITOR CARBON – GRAPHENE OXIDE: N₂ (77 K) AND CO₂ (273 K) ARE COMPLEMENTARY



Zhu, Y.; Murali, S.; Stoller, M.D.; Ganesh, K.J.; Cai, W.; Ferreira, P.J.; Pirkle, A.; Wallace, R.M.; Cychosz, K.A.; Thommes, M.; Su, D.; Stach, E.A.; Ruoff, R.S. *Science* **2011**, 332, 1537-1541

CO₂ MICROPOROUS ANALYSIS OF LI/S BATTERY SUPPORT



- › Microporous carbon support for Li/S batteries
- › Only small micropores (< 1 nm) present in this carbon
- › Therefore, only S₂ (0.4 nm) can be confined within the pores and larger S₄₋₈ molecules (> 0.5 nm) are excluded





PORE SIZE CALCULATIONS

Pore size $\propto p/p_0$

Micropores fill as low as $1 \times 10^{-7} p/p_0$, mesopores $> 0.3 p/p_0$

Traditional: BJH,


Modern: DFT, GCMC






CHEMISORPTION ESSENTIALS

WHAT EXACTLY DOES **CHEMISORPTION** TELL US?



ACTIVE METAL AREA

- › The number of active sites usually expressed as an active metal area




DISPERSION


- › The number of metal atoms accessible to reactive gases expressed as a percentage of the total in the sample

CRYSTALLITE SIZE

- › The approximate size of the metal nanoparticles distributed across the surface of the inert support material




metal atom
support material



CHEMISORPTION ANALYSES

Property	Method
Metal area	Static Isotherms, Flow Pulse Titration
Dispersion & crystallite size	Static Isotherms, Flow Pulse titration
Spillover (excess chemisorption)	Static Isotherms
Oxide reducibility	TPR (Flow)
Catalyst regeneration (carbon removal)	TPO (Flow)
Zeolite acidity	TPD, ammonia (Flow)



ANYTHING ELSE?

CATALYST ACTIVATION

- › Temperature Programmed Reduction: TPR


CATALYST REGENERATION

- › Temperature Programmed Oxidation: TPO

ACID SITE STRENGTH

- › Temperature Programmed Desorption: TPD

ANYTHING ELSE?



CATALYST ACTIVATION


- › Temperature Programmed Reduction: TPR

CATALYST REGENERATION

- › Temperature Programmed Oxidation: TPO

ACID SITE STRENGTH

- › Temperature Programmed Desorption: TPD



ANYTHING ELSE?

CATALYST ACTIVATION

- › Temperature Programmed Reduction: TPR

CATALYST REGENERATION

- › Temperature Programmed Oxidation: TPO

ACID SITE STRENGTH

- › Temperature Programmed Desorption: TPD




HIGH PRESSURE GAS SORPTION ESSENTIALS

ADSORPTION CHARACTERISTICS

› Capture, Storage, Characterization

- **Strong** adsorption
 - Capture
 - Retention
 - Micropores
- **Selectivity**
 - Pore size
 - Adsorption enthalpy





Desirable Properties:

- High adsorption capacity
- Good kinetics
- Recoverable/regenerable

WHY HIGH PRESSURE SORPTION?

Case Study: CO₂ Capture

Activated carbons from anthracite, sub-bituminous coals and lignite



CH₄ STORAGE




Fig. 3 SEM results of selected samples. a WJD, b PY, c DX, d YW, e GQ, f ZC, g ZZ







Fig. 5 Adsorption isotherms of methane on samples

OPTIMIZING HYDROGEN STORAGE FOR CARBON FOOTPRINT REDUCTION



METAL HYDRIDES

PRESSURE-COMPOSITION-TEMPERATURE (PCT PLOT)




The hydrogen moves in and out of the lattice structure, which is completely reversible.

ADSORBENT PERFORMANCE PROPERTIES

HIGH PRESSURE ISOTHERMS & STORAGE CAPACITY

Isotherms On Gravimetric Terms: Adsorption-Desorption
Sample: Activated Carbon Adsorbate: Hydrogen EOS: Helmholtz



BENEFICIAL ADVANCEMENTS

Traditional	Advance	Benefit
Helium void volume	Analysis without helium	Conserve non-renewable resource (helium)
Samples treated as unknown	Intelligent gas dosing based on historical data	High throughput mesopore analysis
Nitrogen/77K micropore analysis	Argon /87K	Time saving ~50%
Liquid argon	87K using LN ₂	Instant availability
25K	Stability at ~20K	Maximizing hydrogen density
Macroscopic models for pore size	Molecular modelling of adsorption behavior in pores of different shapes	Increased pore size accuracy
n/a	Hysteresis scanning	Network information


REVIEW & SUMMARY

- › Surface area: extent of where solids interact with their surroundings
- › Pore size and volume: define how fluids move in, out and through a porous solid
- › Battery electrodes
- › Environmental catalysts
- › Greenhouse gas adsorbents
- › Hydrogen storage
- › All characterized by the various forms of gas sorption measurements



Access Covalent Portals

Customer Access to Data & Community Content



The DATA PORTAL is used by Customers and Lab Partners for uploading and downloading data. It requires two-factor authentication and advanced password protection. Data Portal users have complete access through their home page on the portal to all Community content, and do not require a separate Community account.

[LOGIN TO DATA PORTAL](#)

Community Content Access for **All Other Users**



The COVALENT COMMUNITY PORTAL requires password entry. It contains webinar and other metrology and characterization-related content that we believe would be useful and educational for the materials science innovation community. It does not provide access to any customer data and should only be used by individuals that are not Covalent customers or lab partners.

[LOGIN TO COMMUNITY PORTAL](#)

Or Check Out the New Covalent Academy

55

Covalent Academy

Get more from your data by building your knowledge of advanced materials characterization.

Search resource topics

SEARCH

From
Select

To
Select

Type
Select

Duration
Select

Learning Center

- Browse All Topics

My Learning

Topics In Progress
Completed Topics

LASER ABLATION
INDUCTIVELY
COUPLED PLASMA
MASS
SPECTROSCOPY:
NOT JUST ROCKS

SPEAKER:
Lucas Smith
Director of Business Development
for the Americas,
Teledyne CETAC

February 10, 2022 | 11am PT

TELEDYNE
CETAC TECHNOLOGIES
Laser Microprobe

COVALENT
ACADEMY
Advancements in
Instrumentation Series

Episode 30

Webinar | 60 min

ADVENTURES IN
WAVELENGTH
DISPERSIVE X-RAY
FLUORESCENCE
(WDXRF): FLEXIBLE
ELEMENT ANALYSIS
FOR THIN FILMS
AND MORE

SPEAKER:
Meredith Beebe
Semiconductor X-ray
Metrology Specialist,
Rigaku

January 27, 2022 | 11am PT

Rigaku

COVALENT
ACADEMY
Advancements in
Instrumentation Series

Episode 29

Webinar | 60 min

FAST
CHARACTERIZATION
OF NANOMETER THIN
TO THICK COATINGS
USING PULSED-RF
GLOW DISCHARGE
OPTICAL EMISSION
SPECTROMETRY

SPEAKER:
Philippe Hunault
Technical Sales Elemental
Analysis Specialist,
HORIBA Scientific

December 2, 2021 | 11am PT

HORIBA
Scientific

COVALENT
ACADEMY
Advancements in
Instrumentation Series

Episode 28

Webinar | 60 min

MODERNIZING
MICROSCOPY
METHODS:
CAPABILITIES AND
APPLICATIONS OF
TEM/STEM
SYSTEMS

SPEAKER:
Dr. Jan Ringnald
Principal Scientist, Materials
and Structural Analysis

COVALENT
ACADEMY
Advancements in
Instrumentation Series

Episode 27

Webinar | 60 min

UPGRADING
METROLOGY
SERVICES WITH
MOUNTAINS™ 9:
IMPROVED
AUTOMATION,
VISUALIZATION,
AND ANALYSIS

SPEAKER:
Cyrille Charles
Key Account Manager

COVALENT
ACADEMY
Advancements in
Instrumentation Series

Episode 26

Webinar | 60 min

POROMETRY,
POROSIMETRY, AND
PYCNOMETRY:
THE 3 P'S YOU NEED
FOR POROUS
MATERIALS
CHARACTERIZATION

SPEAKER:
Nanette
Jareenwattananon, PhD
Senior Materials Scientist

COVALENT
ACADEMY
Advancements in
Instrumentation Series

Episode 25

Webinar | 60 min

Watch all episodes
on-demand now at:

<https://academy.covalentmetrology.com>



Q & A Session



**COVALENT
METROLOGY**

Thank You