



**COVALENT
METROLOGY**

Welcome

We'll begin at 11:00 AM PT.

NANOSIMS: HIGH SENSITIVITY IMAGING ANALYSIS FOR DEVICES

SPEAKER:

Junichiro Sameshima, PhD

Senior Manager of
Surface Science Laboratory,
Toray Research Center, Inc.

September 22, 2022 | 11am PT

Toray Research Center, Inc.




COVALENT
ACADEMY

Advancements in
Instrumentation Series

Episode 33



Silicon Valley-based analytical labs and platform delivering quality data and expert analysis for advanced materials and device innovation



Comprehensive Solutions Stack

- 50+ cutting-edge instruments, offering 100+ techniques

- Analytical Services

- Advanced Modeling

- Method Development

- Enterprise Metrology Solutions



Affordable and Fast

- Fast Turnaround Times

- Volume Savings

- Instant Access to Data and Reports in Secure Portal



Flexible Business Model

- Custom Consulting Solutions and Certified On-site Support

- Training and Certification on Instrumentation

- Co-op and Tool-Share Opportunities

- Laboratory Audits



Rich Network of Partnerships

- Partner to the World's Leading Instrument Manufacturers and Labs

- Expanding Instrumentation, Lab Connections, and Learning



Who We Are, Who We Serve

- 500+ Clients, 40 – 60 Added / Quarter

- 60+ Employees, 14 PhDs

- Cutting-edge Analytical Capabilities

- Silicon Valley Lab based in Sunnyvale, CA

PCBA, Semiconductor, and Electronic Device Metrology & Failure Analysis

- DPA / Mechanical Cross-section
- Dye & Pry Test
- EBIC / OBIC failure analysis
- Hot Spot Detection
- IR Imaging / Emission Microscopy
- NIR Imaging
- Root-Cause Failure Analysis

Electron Microscopy and Scanning Probe Microscopy

- AFM & Advanced AFM Modes (EFM, KPFM, MFM, PFM)
- Scanning Acoustic Microscopy (SAM)
- SEM (+ EDS)
- FIB-SEM (+ EDS)
- S/TEM (+ EDS / + EELS)
- Nano-indent / Nano-scratch

Optical Microscopy & Spectroscopy

- Chromatic Aberration
- Digital Optical Microscopy
- FTIR and ATR-FTIR
- Laser Scanning Confocal Microscopy
- Spectral Ellipsometry
- UV-Vis-NIR Spectroscopy
- White Light Interferometry

X-Ray Characterization

- X-Ray Diffraction (XRD)
- X-Ray Reflectometry (XRR)
- Micron-spot ED-XRF
- WDXRF
- Micro-computed X-ray Tomography (Micro-CT)
- 2D X-ray Inspection & X-ray Radiography

Elemental / Chemical Composition Analysis

- EPMA
- GD-OES
- GC-MS
- ICP-MS and LA-ICP-MS
- Raman Microscopy & Spectroscopy
- NMR (1D or 2D; solid / liquid)

Particle Analysis

- Dynamic Light Scattering (DLS)
- Laser Diffraction Particle Size Analysis (PSA)
- Particle Zeta Potential

Material Property Characterization

- DSC
- DMA & TMA
- Rheometry
- TGA
- Surface Zeta Potential
- Porometry / Porosity
- Gas Adsorption
- Gas Pycnometry
- Foam Density
- Tap Density

Surface Spectroscopy Analysis

- Dynamic-SIMS
- ToF-SIMS (Static-SIMS)
- Ion Scattering Spectroscopy (ISS)
- Ultraviolet Photoelectron Spectroscopy (UPS)
- X-ray Photoelectron Spectroscopy (XPS)

Today's webinar is in partnership with

Toray Research Center, Inc.

- Partnership announced in October of 2021
- Toray Research Center, Inc. provides **unique analytical capabilities** for Covalent clients, including:
 - Rutherford Backscattering Spectroscopy (RBS)
 - High-resolution (HR-RBS) and Micro (Micro-RBS)
 - Cathodoluminescence Spectroscopy
 - **Nanoscale Secondary Ion Mass Spectroscopy (NanoSIMS)**

Other Covalent Partners



Dr. Junichiro Sameshima

*Senior Manager of Surface Science Laboratory,
Toray Research Center, Inc.*

- Joined Toray Research Center, Inc. in 1997
- 25 years working with SIMS at Toray – especially for semiconductors research
- Has managed the Surface Science Laboratory since 2012 and obtained the role of Senior Manager in 2018
- Earned Ph.D. for “*Structural analysis of SiC power device and annealing effect by Time-of-Flight Secondary Ion Mass Spectroscopy*” from Osaka University in 2019
 - Master’s degree in Material Physics Engineering from OSAKA University (1997)






High sensitivity imaging analysis for devices using NanoSIMS

Junichiro Sameshima, Ph.D

Surface Science Laboratory,
Toray Research Center, Inc.
Email: junichiro.sameshima.q8.trc.toray

- Introduction to Toray Research Center
- NanoSIMS operational performance
 - Basic performance
 - Advantages over conventional SIMS and other analytical techniques
- Example NanoSIMS applications
 - Semiconductor devices (SiC-MOSFET, 3D NAND flash memory, SRAM)
 - Implanted phosphorous in Si sub. (Quantification)
 - Li-ion Battery (LTO coating, Li-dendrite on carbon negative electrode)
- Wrap-up



Overview

Toray Industries, Inc.

Established in 1926 / Turnover of US\$ is 22.5 billion

Developing and manufacturing advanced materials related to Organic
Synthetic Chemistry, Polymer Chemistry, Nano & Biotechnology

Toray Research Center, Inc.

Established in 1978 / Headquartered in Tokyo, Japan

Labs at Shiga (Central), Kamakura and Nagoya with 530 employees

A professional analytical service provider with 40+ years of accumulated experience

01 Analytical Services
Using latest analytic equipment by
experienced skillful specialists

02 Contract R&D
support

In the fields of Automobile, Displays,
Semiconductor and Advanced
materials., etc.

Locations



UNITED STATES

Fort Wayne, IN
US Agency (Sales)

GERMANY

Munich – AMCEU
(Lab)



CHINA

Shanghai – TRCS (Sales, Lab)



TAIWAN

Taichung – TITP (Sales)

SOUTH KOREA

Anyang – KTS Global (Sales)



JAPAN

Tokyo – Head Office, Sales Office
Shiga – Central laboratory of TRC
Comprehensive Analysis
& Evaluation, Sales Office
Nagoya – Pharmaceutical(GMP)
Sales Office
Kamakura – Bio & Pharmaceutical
(GLP)
Fukuoka – Sales Office

Major Targets

Materials



Semiconductor
& Packaging



Displays
& IT Equipment



Batteries




Automobile
& Aerospace



Life Science

-
- Introduction to Toray Research Center
 - **NanoSIMS operational performance**
 - Basic performance
 - Advantages over conventional SIMS and other analytical techniques
 - Example NanoSIMS applications
 - Semiconductor devices (SiC-MOSFET, 3D NAND flash memory, SRAM)
 - Implanted phosphorous in Si sub. (Quantification)
 - Li-ion Battery (LTO coating, Li-dendrite on carbon negative electrode)
 - Wrap-up

Principles of SIMS




General features of SIMS

Pros:

- All Elements (from H up to U) can be detected (in principle)
- Highest sensitivity (down to ppb levels) among other surface analytical techniques
- Both 2D and 3D ion imaging are available
- Elements contained in the sample can be quantified using a reference Standard
- Isotopic composition analysis is possible due to high mass resolution

Cons:

- Destructive analysis
- “Matrix effect”
Large difference on ionization probability depending on materials and elements



Features of NanoSIMS 50L




- Cosmochemistry
- Geology
- Biology
- Materials
- **Semiconductors / devices**

Specifications

	NanoSIMS 50L	Conventional SIMS	TOF-SIMS
Primary ion	O ⁻ ; for positive ion detection Cs ⁺ ; for negative ion detection	O ₂ ⁺ ; for positive ion detection Cs ⁺ ; for negative ion detection	Bi _n ⁺ * O ₂ ⁺ , Cs ⁺ , GCIB as etching ions
Beam size	50 nm	10 – 100 μm	150 nm – 5 μm
Detection limit	0.1 ppm ~	ppb ~	ppm ~
Mass analyzer	Magnetic Sector	Magnetic Sector or Quadrupole	Time-of-Flight
Number of detected ions	7 ; simultaneously	10 -15	(All ions within mass range)
Depth capability	10 nm – 0.5 μm	1 nm – 50 μm	nm – 10 μm

Capability - Size and Sensitivity -



NanoSIMS 50L


TOF-SIMS

Dynamic SIMS

- XPS : X-ray photoelectron spectroscopy
- AES : Auger electron microscopy
- RBS : Rutherford backscattering spectroscopy
- GD-OES : Glow discharge optical emission spectrometry
- GDMS : Glow discharge mass spectrometry

-
- Introduction to Toray Research Center
 - NanoSIMS operational performance
 - Basic performance
 - Advantages over conventional SIMS and other analytical techniques
 - **Example NanoSIMS applications**
 - Semiconductor devices (SiC-MOSFET, 3D NAND flash memory, SRAM)
 - Implanted phosphorous in Si sub. (Quantification)
 - Li-ion Battery (LTO coating, Li-dendrite on carbon negative electrode)
 - Wrap-up

Commercially available SiC-MOSFET*
Silicon carbide metal-oxide-semiconductor field-effect transistor




Cross-section was prepared by FIB process

NanoSIMS Imaging for SiC-MOSFET cross-section

Positive ions


Sputtered by O⁻ primary beam



- Al: p-well
- B: poly-Si electrode

Negative ions


Sputtered by Cs⁺ primary beam




- P: PSG and n-source

Comparison between NanoSIMS and TEM-EDX imaging for dopants


HAADF-STEM




NanoSIMS




Almost no dopant was detected at all (except for P in PSG)






B



Al




P




Dopant distributions are clearly captured

Al distribution in p-well region observed by STEM and NanoSIMS

BF-STEM Image



NanoSIMS Image




- Aluminum ion distribution by NanoSIMS corresponds to crystal defect generated by ion implantation.

Comparison between NanoSIMS and SCM

SCM Image


SCM : Scanning Capacitance Microscopy

- Active carrier distribution
- p-n junction




NanoSIMS Image

- Dopant species
- Dopant distribution




Stack structure for NAND flash memory

Commercially available 3D memory




Negative ions detected
Sputtered by Cs⁺ primary beam

Backside NanoSIMS Image




- F detected in W (WF) layer - - > originated from precursor in CVD process

Quantification and verification of detection limit - P in Si wafer -



Quantified line profile of Phosphorous in Si

NanoSIMS line profiling



For Si device sample within less than 2 μm area

———— FIB + NanoSIMS —————>


Profiles @ detection limit of 10¹⁵ atoms/cm³ levels

-
- Introduction to Toray Research Center
 - NanoSIMS operational performance
 - Basic performance
 - Advantages over conventional SIMS and other analytical techniques
 - Example NanoSIMS applications
 - Semiconductor devices (SiC-MOSFET, 3D NAND flash memory, SRAM)
 - Implanted phosphorous in Si sub. (Quantification)
 - Li-ion Battery (LTO coating, Li-dendrite on carbon negative electrode)
 - Wrap-up




Active material for Li-ion battery: evaluation of LTO coating

Positive ions
Sputtered by O⁻ primary beam

Imaging analysis of LTO layer (Field of view: 10 μm)






LTO thick



Uniformity of LTO coating was evaluated based on Ti and Co distributions.

LTO thin



Low intensity of Ti

Secondary ion image of Li, Ti, Co

“LTO thick” - Almost all LCO particles are covered with the LTO film.





“LTO thin” - Some LCO particles are not covered with the LTO film or the film thickness is thinner than the surrounding particles.

Li-dendrite measurement on carbon negative electrode – Li-ion Battery -

Experimental flow

- Battery disassembly in an inert atmosphere
↓
- Cross-section processing in an inert atmosphere (Cross-section polishing)
↓
- Transport to instrument in an inert environment
↓
- Mapping by NanoSIMS

FIB-SEM (Cross-Sectional view)




Evaluation of Li, which is hard to detect by EDX, is enabled by NanoSIMS


Comparison – NanoSIMS, SEM-EDX, EPMA -

STD: Standard Materials

	NanoSIMS 50L	SEM-EDX	EPMA
Probe	Ion beam (Cs ⁺ , O ⁻)	Electron beam	Electron beam
Special resolution	~50 nm	20~100 nm	~100 nm
Detection limit	~ppm	sub wt%	0.05~0.1 %
Detectable elements	H~U	B~U	B~U
Qualitative analysis	No (7 elements have to decided in advance)	Yes	Yes
quantification	Yes (STD needed)	Yes (for semi-quantification, STD not needed)	
Other Pros.	- Highest sensitivity	- short measurement time	- Higher sensitivity than SEM-EDX



NanoSIMS is capable of detecting Si along with grain boundary in ceramics, contrary to SEM-EDX, EPMA



-
- Introduction of Toray Research Center
 - NanoSIMS performance
 - Basic performance
 - Advantage towards conventional SIMS or other analytical techniques
 - NanoSIMS applications
 - Semiconductor devices (SiC-MOSFET, 3D NAND flash memory, SRAM)
 - Implanted phosphorous in Si sub. (quantification)
 - Li-ion Battery (LTO coating, Li-dendrite on carbon negative electrode)
 - Organic Coating layer
 - **Wrap-up**

Summary

- NanoSIMS is optimized for **high sensitivity elemental imaging** analysis.
- NanoSIMS can provide
 - **2D and 3D Imaging**
 - Line profiles (which can be quantified with a reference sample)
- NanoSIMS can *measure* dopant / impurities contained at *lower levels*, which are invisible to TEM-EDX, SEM-EDX, EPMA
- NanoSIMS often works better by combining it with other analytical methods.

Summary

- NanoSIMS is a great choice in case you need:
 - dopant imaging [implantation, diffusion, additive, etc.]
 - impurities imaging [contamination, segregation, coexistence, etc.]
 - depth profiles [p-/n-type dopant in semiconductor]on **devices, small target areas, powders**, etc.

- Toray Research Center has developed not only the NanoSIMS analytical technique, but also sample preparation for analyses across manifold fields for various analytical methods. By taking advantage of this, we can get out of the best NanoSIMS performance for the field below
 - **Semiconductor devices**
 - **Li-ion Battery, SOFC**
 - **Materials including layer films, powder**
 - **Tissues, Cells [not shown in this seminar]**

Toray Research Center hopes to help promote your R&D with you.



Thank you for your attention.

Stay Tuned!

We'll announce the next episode soon on our website at:

<https://covalentmetrology.com>

Thank you for attending!

To show our appreciation, we're offering all attendees a **special limited-time discount**

Up to 15% Attendee Discount*

On Your Next NanoSIMS Project

**10% discount for jobs quoted at \$1,500 - \$2,999; jobs over \$3,000 qualify for the 15% discount*

Schedule a **Free Discovery Consult Appointment** to discuss your needs with a Covalent expert

[Calendly Link in the Chat](#)

All Episodes On-Demand in the Covalent Academy

Covalent Academy

Get more from your data by building your knowledge of advanced materials characterization.

Search resource topics **SEARCH**

From Select To Select Type Select Duration Select

Learning Center

- [Browse All Topics](#)

My Learning

Topics In Progress

Completed Topics

Find all our past episodes at

academy.covalentmetrology.com

LASER ABLATION INDUCTIVELY COUPLED PLASMA MASS SPECTROSCOPY: NOT JUST ROCKS

SPEAKER:
Lucas Smith
Director of Business Development for the Americas, Teledyne CESC

February 10, 2022 | 1 Item PT

COVALENT METROLOGY

COVALENT ACADEMY
Advancements in Instrumentation Series
Episode 30

Webinar | 60 min

ADVENTURES IN WAVELENGTH DISPERSIVE X-RAY FLUORESCENCE (WDXRF): FLEXIBLE ELEMENT ANALYSIS FOR THIN FILMS AND MORE

SPEAKER:
Meredith Beebe
Semiconductor X-ray Metrology Specialist, Rigaku

January 27, 2022 | 1 Item PT

COVALENT METROLOGY

COVALENT ACADEMY
Advancements in Instrumentation Series
Episode 29

Webinar | 60 min

FAST CHARACTERIZATION OF NANOMETER THIN TO THICK COATINGS USING PULSED-RF GLOW DISCHARGE OPTICAL EMISSION SPECTROMETRY

SPEAKER:
Philippe Hunault
Technical Sales Elemental Analysis Specialist, HORIBA Scientific

December 2, 2021 | 1 Item PT

COVALENT METROLOGY

COVALENT ACADEMY
Advancements in Instrumentation Series
Episode 28

Webinar | 60 min

MODERNIZING MICROSCOPY METHODS: CAPABILITIES AND APPLICATIONS OF TEM/STEM SYSTEMS

SPEAKER:
Dr. Jan Ringsdorf
Principal Scientist, Materials and Structural Analysis, Thermo Fisher Scientific

November 9, 2021 | 1 Item PT

COVALENT METROLOGY

COVALENT ACADEMY

UPGRADING METROLOGY SERVICES WITH MOUNTAINS™ 9: IMPROVED AUTOMATION, VISUALIZATION, AND ANALYSIS

SPEAKER:
Cyrille Charles
Key Accounts Manager, Digital Surf

October 21, 2021 | 1 Item PT

COVALENT METROLOGY

POROMETRY, POROSIMETRY, AND PYCNOMETRY: THE 3 P'S YOU NEED FOR POROUS MATERIALS CHARACTERIZATION

SPEAKER:
Nanette Jarenwattananon, PhD
Senior Manager, Material Property Testing

October 7, 2021 | 1 Item PT

COVALENT METROLOGY



Q & A Session



**COVALENT
METROLOGY**

Thank You

For Q and A

Sample Preparation for cross-section for semiconductor or materials

The way to make sample	Sample	Pros.	Concerns and solution
FIB preparation	Specific area in device	Pick up an interested are	Damage due to Ga ion beam -- > pre-sputtering
Broad Ion Beam Etching	Materials	Clean	Curved sample -- > keep wide area
cleavage	Semiconductor wafer	As it is, intact	Roughness or line at surface -- > additional polish etc.

NanoSIMS 50L synopsis

