



**COVALENT
METROLOGY**

Welcome

covalentmetrology.com

SAMPLE PREPARATION ORIENTATION AND ITS VALUE FOR TRANSMISSION ELECTRON MICROSCOPY (TEM)

SPEAKER:


Ryan Dudschus

Senior Metrology Engineer,
Covalent Metrology

July 15, 2021 | 11AM PT

COVALENT
ACADEMY


Episode 24





COVALENT METROLOGY

Silicon Valley-based analytical labs and platform delivering quality data and expert analysis for advanced materials and device innovation



Comprehensive Solutions Stack

*50+ cutting-edge
instruments, offering
100+ Techniques*

Analytical Services

Advanced Modeling

Method Development

Temp. Staffing Solutions



Affordable and Fast

*Fast Turnaround Times,
No Expedite Fees*

Volume Savings

*Instant Access to Data
and Reports in Secure
Portal*



Flexible Business Model

*Custom Consulting
Solutions and Certified
Onsite Support*

*Training and Certification
on Instrumentation*

*Co-op and Tool-Share
Opportunities*


Laboratory Audits



Rich Network of Partnerships

*Partner to World's
Leading Instrument
Manufacturers and Labs*

*Expanding
Instrumentation, Lab
Connections and
Learning*



Who We Are, Who We Serve

*500 Clients + 40-60 New
Clients/ Quarter*

40 People, 13 PhDs

*Cutting-edge Analytical
Capabilities*

*Lab Locations:
Sunnyvale, CA*

Covalent Analytical Service Categories

4

PCBA, Semiconductor, and Electronic Device Metrology & Failure Analysis

- DPA / Mechanical Cross-section
- Dye & Pry Test
- EBIC / OBIC failure analysis
- Hot Spot Detection
- IR Imaging / Emission Microscopy
- NIR Imaging
- Root-Cause Failure Analysis

Electron Microscopy and Scanning Probe Microscopy

- AFM & Advanced AFM Modes (EFM, KPFM, MFM, PFM)
- Scanning Acoustic Microscopy (SAM)
- SEM (+ EDS)
- FIB-SEM (+ EDS)
- S/TEM (+ EDS / + EELS)
- **Nano-indent / Nano-scratch**

Optical Microscopy & Spectroscopy

- Chromatic Aberration
- Digital Optical Microscopy
- **FTIR and ATR-FTIR**
- Laser Scanning Confocal Microscopy
- Spectral Ellipsometry
- UV-Vis-NIR Spectroscopy
- White Light Interferometry

X-Ray Characterization

- X-Ray Diffraction (XRD)
- X-Ray Reflectometry (XRR)
- Micron-spot ED-XRF
- WDXRF
- Micro-computed X-ray Tomography (Micro-CT)
- 2D X-ray Inspection & X-ray Radiography

Elemental / Chemical Composition Analysis

- EPMA
- **GD-OES**
- GC-MS
- **ICP-MS and LA-ICP-MS**
- **Raman Microscopy & Spectroscopy**
- NMR (1D or 2D; solid / liquid)

Particle Analysis

- Dynamic Light Scattering (DLS)
- Laser Diffraction Particle Size Analysis (PSA)
- Particle Zeta Potential

Material Property Characterization

- DSC
- DMA & TMA
- Rheometry
- TGA
- Surface Zeta Potential

Coming Soon:

- *Porometry / Porosity*
- *Gas Adsorption*
- *Gas Pycnometry*
- *Foam Density*
- *Tap Density*

Surface Spectroscopy Analysis

- Dynamic-SIMS
- ToF-SIMS (Static-SIMS)
- Ion Scattering Spectroscopy (ISS)
- Ultraviolet Photoelectron Spectroscopy (UPS)
- X-ray Photoelectron Spectroscopy (XPS)

Ryan Dudschus, PhD



Senior Metrology Engineer, Lamella Preparation
Covalent Metrology

- Senior metrology engineer working primarily in Electron and Scanning Probe Microscopy groups
 - Using mainly SEM and AFM techniques
- He earned his B.S. and M.S. in Physics from UC Santa Cruz
 - Research focused in condensed matter
- Ryan's M.S. research used the Stanford Synchrotron Radiation Lightsource (SSRL) for Extended X-ray Absorption Fine Structure (EXAFS) data collection (instead of an SEM, XRD apparatus, or other standard measurement device)



Sample Prep. Correlation to Quality of TEM Image	7
Fundamentals	8
Grid Preparation and Types of Lamella Thinning	15
Sample Orientation and Variety of Lamella Lift-out Procedures	21
Lamella Defects	30
Conclusion	39

Sample preparation accounts for ~90% of the TEM image quality...



Dualbeam Fundamentals for TEM Sample Preparation


Basic DualBeam setup

SEM use cases:

- Imaging
- Depositions


FIB use cases:

- Milling/Thinning
- Deposition




Lamella Preparation Overview – Step by Step

10




Thinner lamella provides in higher quality TEM imaging




Physical example of TEM imaging between 100nm thick and 40nm thick lamella

~100 nm thick




~100 nm
thick

~40 nm
thick




~40 nm thick




Note: the thinner the lamella, the more challenging further thinning becomes



- Ga FIB has an interaction volume based on its lining energy.
- Ga ions may knock atoms off of lattice site.
- Ga ions may imbed themselves in the specimen



[2]




- Ga FIB has **gaussian shape**
- Each material has own mill rate
- Over-tilt affects the shape of the lamella
- Lower currents for precise mills
- Higher currents for faster mills
- Short dwell time for uniform, accurate pattern shapes
- High dwell time for mill harder materials



Grid Preparation and Types of Lamella Thinning

- Grid are chosen based on material composition of the specimen
- Grids could be made from various materials:
Au, Be, C, Cu, Si, Mo, Ni

Ni – Mesh Grid



Narrow Posts



Cu




Wide Posts



Mo




Grid Preparation



- Creates a clean surface for lamella attachment.
- Top mount “T-pit” prevents redeposition on lamella from grid.

Lamella Attachment Points




- **Top Mount**
 - Two points of contact. Less prone to twisting (but curling still possible)
- **Side Mount**
 - Quick and more secure attachment; speeds up time to make a lamella
 - More prone to bending/twisting

Basic Types of Thinning Lamella


18

Windowed




- Thick bottom helps prevent lamella bending.
- Lamella curling may still occur.

Pillar



- Good for near-surface layers/features.
- Pillar provides stabilization.
- Lamella bending and curling may still occur.

Through






- Useful for large multi-layers/features.




Lamella – Electron Transparency

19




10kV



5kV




2kV



Decreasing Thickness →

Note: electron transparency thickness is not the same for all materials

10kV Transparency




2kV Transparency


Lamella – Electron Transparency Thickness

20



10kV



5kV




2kV



Sample Orientation and Variety of Lamella Lift-out







Sample Orientation Planning

22




Lamella Lift-out: Top-Down – Process

23




Lamella Lift-out: Top-Down

24











- Good for near-surface features
- Smooth top surfaces are ideal
- Most common prep technique




Lamella Lift-out: Inverted – Process

25




Lamella Lift-out: Inverted

26



Used for:


- Samples with a rigid/rough top surface
- Substantially harder material near lamella top to use lower overtilt angle and prevent milling through softer material before harder material is substantially thinner.





Lamella Lift-out: Planar – Process

27


1




2




3



4







5



Lamella Lift-out: Planar


28



Variety of Lamella Lift-out – Recap


29

Top-Down




- Most common lamella type

Inverted



- Provide a uniform layer above feature to prevent curtaining
- Place lower mill rate material near top to use a lower overtilt angle while thinning.

Planar



- Examines layer/features on the x-y plane.
- Isolate individual layers for TEM/EDS analysis
- Provide another view on devices (planar cross-section)

Lamella Defects

Lamella Defects – Amorphous Damage




31

Amorphous damage in Si due to:

1a) 30 kV, ~ 21.5 nm

1b) 5 kV, ~ 6.6 nm

1c) 2 kV, ~ 3.1 nm






TEM imaged of Si polished with:

2a) 30 kV

2b) 5 kV


2c) 2 kV




[3]


Thank you to Mark Najarian at ThermoFisher

Reduction in TEM image quality due to variation in thickness





Redeposition





- Lamella may bend at various thicknesses based on the sample material





Thinned lamella region
bending and curling
upwards

- Deformations in feature/boundary interface can lead to inaccurate measurements

Lamella Defects – Curtaining


34





Curtaining occurs in:

- **Soft materials**
- Boundaries with a change in texture
- Grains
- **Devices/features**

- Incorrect over-tilt angles – or specimens with materials that have different milling rates – result in “wedged” lamella
- Wedge lamella have the thickness either increase or decrease across the thinned window
- Results in loss of TEM image resolution across the thinned window








Holes may contribute to other lamella defects, such as bending/curling, as well as possibly destroying the feature of interest.





Lamella Defects – TEM Images

37



Defects Don't Have to Stop a Good TEM Image

38



- Lamella preparation is not only science, but a form of art
- Sample preparation is crucial to TEM imaging
- Understand the material being processed
- Select grids that do not contain materials of the specimen
- Sample orientation
- Several different way to thin a lamella
- Electron transparency differs amongst lining energy and material type
- Lamella defects have numerous root causes negatively impacting TEM image quality

Future webinars about TEM sample preparation will focus on advanced lamella thinning techniques such as cross-section lamella for TEM imaging, rocking mill, targeted samples, etc.

1. Brüggemann, D., Wolfrum, B., & de Silva, J. P. (2014). Fabrication, properties and applications of gold nanopillars. In *Handbook of nanomaterials properties* (pp. 317-354). Springer, Berlin, Heidelberg.
 - FIB-SEM diagram
2. Shahali, Hesam & Hasan, Jafar & Wang, Hongxia & Tesfamichael, Tuquabo & Yan, Cheng & Yarlagadda, Prasad. (2019). Evaluation of Particle Beam Lithography for Fabrication of Metallic Nano-structures. *Procedia Manufacturing*. 30. 261-267. 10.1016/j.promfg.2019.02.038.
 - FIB interaction volume
3. Giannuzzi, L. A., Geurts, R., & Ringnalda, J. (2005). 2 keV Ga⁺ FIB milling for reducing amorphous damage in silicon. *Microscopy and Microanalysis*, 11(S02), 828-829.
 - Amorphous damage layer thickness
 - Provided by Mark Najarian @ ThermoFisher
4. “Gaussian Beam Propagation.” *Edmund Optics*, www.edmundoptics.eu/knowledge-center/application-notes/lasers/gaussian-beam-propagation/.
 - Gaussian beam image

- Thank you to **Roozbeh Nikkah-Moshaie** and **Yi Zhang** for provided the TEM images
- Thank you to **Aleia Kim** for the graphics and formatting
- Thank you to **Covalent Metrology Services** for making this webinar possible

Stay Tuned! The next episode will be announced soon at our Website:

www.covalentmetrology.com

And on LinkedIn: www.linkedin.com/company/covalentmetrology

Want to learn more about Covalent's
Atomic Force Microscopy services?

Talk with a Covalent Expert!

[Schedule your Appointment Now with Calendly](#)

- link is in the chat -



Science Forward.

Covalent delivers quality data and expert analysis for advanced materials and device innovation.

Type technique or measurement to find metrology solutions

[All Analytical Services](#) ▾




 [Learn more about Advanced Search Beta](#)

[Contact Us](#)

Access Covalent Portals

Customer Access to Data & Community Content



The DATA PORTAL is used by Customers and Lab Partners for uploading and downloading data. It requires two-factor authentication and advanced password protection. Data Portal users have complete access through their home page on the portal to all Community content, and do not require a separate Community account.

LOGIN TO DATA PORTAL

Community Content Access for **All Other Users**



The COVALENT COMMUNITY PORTAL requires password entry. It contains webinar and other metrology and characterization-related content that we believe would be useful and educational for the materials science innovation community. It does not provide access to any customer data and should only be used by individuals that are not Covalent customers or lab partners.

LOGIN TO COMMUNITY PORTAL



Q & A Session



**COVALENT
METROLOGY**

Thank you

covalentmetrology.com